



**OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT
FACILITIES DEVELOPMENT DIVISION**

**APPLICATION FOR OSHPD PREAPPROVAL
OF MANUFACTURER'S CERTIFICATION (OPM)**

OFFICE USE ONLY

APPLICATION #: OPM-0067-13

OSHPD Preapproval of Manufacturer's Certification (OPM)

Type: New Renewal Update to Pre-CBC 2013 OPA Number: _____

Manufacturer Information

Manufacturer: Entermics Medical Systems

Manufacturer's Technical Representative: Adam Van Essen

Mailing Address: W164 N9283 Water Street, Menomonee Falls, WI 53052

Telephone: 262-509-6457

Email: adame@entermics.com

Product Information

Product Name: EC Series Warming Cabinets

Product Type: Warming Cabinets

Product Model Number: EC770, EC770L, EC1260BL, EC1540, EC1540BL, EC1730BL, and EC2060

General Description: Warming cabinets restraint attachment assemblies 5017207 and 5017226 secured to various adjacent support structures.

Applicant Information

Applicant Company Name: Entermics Medical Systems

Contact Person: Adam Van Essen

Mailing Address: W164 N9283 Water Street, Menomonee Falls, WI 53052

Telephone: 262-509-6457

Email: adame@entermics.com

I hereby agree to reimburse the Office of Statewide Health Planning and Development review fees in accordance with the California Administrative Code, 2013.

Signature of Applicant: *Adam Van Essen*

Date: 12/3/2013

Title: Quality Assurance Manager

Company Name: Entermics Medical Systems

"Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"

STATE OF CALIFORNIA – HEALTH AND HUMAN SERVICES AGENCY
OSH-FD-700 (REV 5/30/13)





**OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT
FACILITIES DEVELOPMENT DIVISION**

Registered Design Professional Preparing Engineering Recommendations

Company

Name: IDS Group, Inc.

Name: Rami Elhassan California License Number: SE3930

Mailing Address: 1 Peters Canyon Rd., Suite 130, Irvine, CA 92606

Telephone: 949-387-8500 Email: Rami.Elhassan@idsgi.com

OSHPD Special Seismic Certification Preapproval (OSP)

- Special Seismic Certification is preapproved under OSP- (Separate application for OSP is required)
- Special Seismic Certification is no preapproved

Certification Method(s)

- Testing in accordance with: ICC-ES AC156 FM 1950-10
- Other* (Please Specify): _____


*Use of criteria other than those adopted by the California Building Standards Code, 2013 (CBSC 2013) for component supports and attachments are not permitted. For distribution system, interior partition wall, and suspended ceiling seismic bracings, test criteria other than those adopted in the CBSC 2013 may be used when approved by OSHPD prior to testing.

- Analysis
- Experience Data
- Combination of Testing, Analysis, and/or Experience Data (Please Specify): _____

List of Attachments Supporting the Manufacturer's Certification

- Test Report Drawings Calculations Manufacturer's Catalog
- Other(s) (Please Specify): Manufacturer's specification sheets

OFFICE USE ONLY – OSHPD APPROVAL VALID FOR CBC 2013 ONLY

Signature:  Date: June 06, 2014

Print Name: Jeffrey Y. Kikumoto

Title: Senior Structural Engineer

Condition of Approval (if applicable): _____

"Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"



Office of Statewide Health, Planning and Development
 Pre-Approval of Manufacture's Certification (OPM)

OPM-0067-13

EQUIPMENT MANUFACTURER: Entermics Medical Systems

EQUIPMENT TYPE: Warming Cabinets

Cabinet Models: 7-13

EC 770

EC 770 L

EC 1260 BL


EC 1540

EC 1540 BL

EC 1730 BL

EC 2060

EC 2140 BL



SHEET TITLE:
COVER SHEET

EQUIPMENT INFORMATION:
 MODELS EC770, EC770L, EC1260BL, EC1540, EC1540BL, EC1730BL, EC2060, EC2140 BL




ISSUE DATE:
 JUNE 6, 2014

OPM SHEET NO:
0.1

ENG STAMP:

OPM INDEX	
SHEET NUMBER	SHEET NAME
0.1	COVER SHEET
0.2	DRAWING INDEX/ FLOWCHART
0.3	NARRATIVE AND DESIGN PROFESSIONAL NOTES
1.1	GENERAL NOTES
1.2	GENERAL NOTES
1.3	GENERAL NOTES
1.4	GENERAL NOTES
1.5	GENERAL NOTES
2.1	CABINET CONTENTS AND WEIGHTS
2.2	CABINET CONFIGURATION AND DIMENSIONS
2.3	NUMBER OF ATTACHMENT BRACKETS REQUIRED
2.4	4 ATTACHMENTS - CABINET TOP VIEW AND SIDE ELEVATION
2.5	6 OR 8 ATTACHMENTS - CABINET TOP VIEW AND SIDE VIEW
2.6	REFERENCE BRACKET ASSEMBLY TYPE A
2.7	REFERENCE BRACKET ASSEMBLY TYPE B
3.1	ATTACHMENT TO METAL STUD WALL
3.2	BRACKET TO SHEET METAL
3.3	BACKING TRACK CONFIGURATION
4.1	ATTACHMENT TO CONCRETE WALL
4.2	BRACKET TO CEMENTITIOUS MATERIAL
5.1	ATTACHMENT TO CMU WALL
6.1	BRACKET TO SURFACE MOUNTED PLATE
6.2	SURFACE PLATE CONFIGURATION
7.1	ATTACHMENT TO WOOD STUD WALL
7.2	BRACKET TO WOOD
7.3	WOOD BACKING CONFIGURATION
8.1	BRACKET TYPE A OR B ATTACHMENT TO BACK OF CABINET
8.2	BRACKET TYPE A OR B ATTACHMENT TO SIDE OF CABINET

DATE: 06/06/2014



SHEET TITLE:
DRAWING INDEX/ FLOWCHART

EQUIPMENT INFORMATION:
MODELS EC770, EC770L, EC1260BL, EC1540, EC1540BL, EC1730BL, EC2060, EC2140 BL



ISSUE DATE:
JUNE 6, 2014

OPM SHEET NO:
0.2

ENG STAMP:

NARRATIVE

THE FOLLOWING OPM DOCUMENT DEFINES THE REQUIRED STRUCTURAL RESTRAINT FOR THE NOTED WARMING CABINETS. PLEASE SEE THE FOLLOWING GENERAL NOTES THAT OUTLINE THE DESIGN CRITERIA FOR THIS OPM DOCUMENT.

A FABRICATED ATTACHMENT BRACKET IS UTILIZED NEAR EACH OF THE FOUR WARMING CABINET CORNERS FOR ATTACHMENT TO ADJACENT SUPPORT STRUCTURE TO PROVIDE LATERAL RESTRAINT FOR THE CABINET. THE CUSTOM-FABRICATED ATTACHMENT ALLOWS FOR MOUNTING LOCATION FLEXIBILITY TO BOTH THE CABINET AND SUPPORT STRUCTURE.

THE BRACKETS AND THE HARDWARE NEEDED TO ATTACH THE BRACKET TO BOTH THE CABINET AND TO THE EXTENSION (FOR TYPE B ONLY SHOWN ON SHEET 2.7) AS REQUIRED BY THIS DOCUMENT IS PROVIDED BY ENTHERMICS. THE SCREWS AND ANCHORS NEEDED TO CONNECT THE BRACKETS TO THE WALL AS REQUIRED BY THIS DOCUMENT SHALL BE PROVIDED BY THE CONTRACTOR. THE INSTALLATION OF THE CABINETS SHALL BE BY THE CONTRACTOR.

RESPONSIBILITIES OF THE STRUCTURAL ENGINEER OF RECORD

THE STRUCTURAL ENGINEER OF RECORD (SEOR) SHALL:

1. VERIFY THE ADEQUACY OF THE STRUCTURE (SUCH AS FLOORS, BEAMS, WALLS, ETC.) WHICH SUPPORT/RESTRAINT THE UNIT FOR THE LOADS/FORCES IMPOSED ON THEM BY THE UNIT AS WELL AS ALL OTHER LOADS.
2. PROVIDE ANY SUPPORTING STRUCTURE REQUIRED TO SUPPORT/RESTRAINT THE WEIGHTS AND FORCES SHOWN, IN ADDITION TO ALL OTHER LOADS.
3. VERIFY THAT THE INSTALLATION IS IN CONFORMANCE WITH THE 2013 CALIFORNIA BUILDING CODE AND WITH THE DETAILS SHOWN IN THIS PRE-APPROVAL. VERIFY THAT THE ACTUAL EQUIPMENT'S WEIGHT, CG LOCATION, ANCHOR LOCATIONS, ANCHOR DETAILS AND THE MATERIAL AND GAGE OF THE UNIT WHERE ATTACHMENTS ARE MADE AGREE WITH THE INFORMATION SHOWN ON THE PRE-APPROVAL DOCUMENTS.
4. SPECIFY THIS PRE-APPROVAL:
 - a. IDENTIFY THE CABINETS BASED UPON THE JOB SPECIFIC CRITERIA.
 - b. IDENTIFY THE ADJACENT SUPPORT STRUCTURE FOR EACH CABINET AND METHOD OF ATTACHMENTS.
 - c. IDENTIFY THE CORRECT RESTRAINT ATTACHMENT DETAIL FOR EACH CABINET BASED UPON THE JOB SPECIFIC SEISMIC CRITERIA, AND LOCATION WITHIN BUILDING.
5. PROVIDE CONSTRUCTION DOCUMENTS AS NEEDED TO PROPERLY CONVEY THE MATERIALS AND DETAILS NEEDED FOR PROPER SYSTEM INSTALLATION. CONSTRUCTION DOCUMENTS SHALL INCLUDE THE FOLLOWING AT A MINIMUM:
 - a. PLAN SHOWING CABINET LOCATION WITHIN BUILDING.
 - b. SITE SPECIFIC DESIGN CRITERIA.
 - c. PLAN VIEW OF ADJACENT STRUCTURE PROVIDING CABINET RESTRAINT.
 - d. LOCATION OF ANY SUPPLEMENTAL STRUCTURE SUPPORT/STRENGTHENING.
 - e. OPM SHEET REFERENCE FOR ALL DETAILS USED.

 ENG STAMP:	SHEET TITLE: NARRATIVE AND DESIGN PROFESSIONAL NOTES	
	EQUIPMENT INFORMATION: MODELS EC770, EC770L, EC1260BL, EC1540, EC1540BL, EC1730BL, EC2060, EC2140 BL	
	 1 Peters Canyon Rd. Irvine, CA 92606 www.idsgi.com 949.387.8500	

GENERAL NOTES

1. THIS OSHPD PREAPPROVAL OF MANUFACTURER'S CERTIFICATION (OPM) IS BASED ON THE CBC 2013. THE DEMAND (DESIGN FORCES) FOR USE WITH THIS OPM SHALL BE BASED ON THE CBC 2013.

2. THE CONTRACTOR SHALL NOTIFY OSHPD, INSPECTOR OF RECORD (IOR), AND THE STRUCTURAL ENGINEER OR RECORD (SEOR) WHERE A CONFLICT OR DISCREPANCY OCCURS BETWEEN THE OPM AND ANY OTHER PORTION OF THE CONSTRUCTION DOCUMENTS, FIELD CONDITIONS OR WHERE ANY CONDITIONS ARISE NOT COVERED BY THESE DOCUMENTS WHEREIN WORK WILL NOT COMPLY WITH THE REQUIREMENTS OF THE CODES NOTED IN NOTE 1 ABOVE.

3. DESIGN CRITERIA:

a. CODE: 2013 CALIFORNIA BUILDING CODE (CBC).

b. LOAD COMBINATIONS: ALLOWABLE STRENGTH DESIGN (ASD) IN ACCORDANCE WITH ASCE 7-10 SECTION 12.4

c. SYSTEM WEIGHTS:

CABINET SELF WEIGHT FROM MANUFACTURER'S LITERATURE

BLANKET WEIGHT IS 20 POUNDS PER CUBIC FOOT OF BLANKET AREA VOLUME IN THE CABINET PER 2013 CBC SECTION 1607A, TABLE 1607A.1 FOR OTHER MEDIA.

IV FLUID WEIGHT IS 2.2 POUNDS PER LITER

d. USE OF ITEMS WITH REFERENCED RESEARCH REPORTS IS BASED UPON EACH REPORT'S SECTION "5.0 CONDITIONS OF USE"

4. SEISMIC CRITERIA:

a. THE DETAILS IN THIS PRE-APPROVAL MAY BE USED AT ANY LOCATION IN THE STATE OF CALIFORNIA WHERE $Sds \leq 2.5$.

b. FORCES PER ASCE 7-10 SECTION 13.3.1 (EQUATIONS 13.3-1, 13.3-2, 13.3-3)

WHERE:

$a_p = 1.0$ (ASCE 7-10, TABLE 13.5-1, CABINETS).

$R_p = 2.5$ (ASCE 7-10, TABLE 13.5-1, CABINETS).

$I_p = 1.5$ (CBC 2013 SECTION 1616A.1.17).

$Z/H = \text{VARIABLE} \leq 1.0$ (DETERMINED BY INSTALLATION WITHIN BUILDING).


$\rho = 1.0$ (ASCE 7-10, SECTION 13.3.1).

$\Omega_0 = 1.0$ (EXCEPT FOR ANCHORAGE TO CONCRETE = 2.5 PER ASCE 7-10 TABLE 13.5-1, CABINETS), (ASCE 7-10, SECTION 13.3.1).

$Sds = \text{VARIABLE} \leq 2.5$ (DETERMINED BY SPECIFIC LOCATION SEISMICITY).

$F_v = 0.2 * Sds * D$ (ASCE 7-10, SECTION 13.3.1), APPLIED CONCURRENTLY WITH F_p .

$F_p = \text{SEISMIC DESIGN FORCE IN ACCORDANCE WITH ASCE 7-10 SECTION 13.3.1. THIS PRE-APPROVAL INCORPORATES } F_p \text{ VALUES BASED ON THE SITE SPECIFIC } Sds \text{ AND } Z/H \text{ VARIABLES NOTED ABOVE. SEE SHEET 2.3 FOR THE NUMBER AND TYPE OF ATTACHMENT BRACKETS REQUIRED.}$



ENG STAMP:

SHEET TITLE:

GENERAL NOTES

EQUIPMENT INFORMATION:

MODELS EC770, EC770L, EC1260BL, EC1540, EC1540BL, EC1730BL, EC2060, EC2140 BL



ISSUE DATE:

JUNE 6, 2014

OPM SHEET NO:



1.1

5. STRUCTURAL STEEL:

- a. STEEL MATERIALS SHALL CONFORM TO THE FOLLOWING, UNLESS OTHERWISE NOTED (U.N.O.):
 PLATES AND BARS: ASTM A36 (Fy=36ksi, Fu=58ksi)
 BOLTS: ASTM A307, U.N.O.
 NUTS: ASTM A563
 WASHERS: ASTM F436
 STAINLESS STEEL: ASTM A240 TYPE 304 OR BETTER (Fy=30ksi, Fu=75ksi)
- b. ALL STEEL (EXCEPT STAINLESS STEEL) SHALL BE SHOP PRIMED WITH ZINC OXIDE PRIMER UNLESS OTHERWISE NOTED.
- c. PAINT ALL STRUCTURAL STEEL WITH WEATHER/RUST RESISTANT PAINT IN ACCORDANCE WITH THE PROJECT SPECIFICATIONS AND ARCHITECTURAL DETAILS UNLESS OTHERWISE NOTED.
- d. DIAMETER OF BOLT HOLE SHALL BE 1/16" LARGER THAN THE BOLT'S DIAMETER UNLESS OTHERWISE NOTED.
- e. PROVIDE A STANDARD WASHER AT ALL TURNING SURFACES. LOCK WASHERS MAY BE USED IN LIEU OF STANDARD WASHERS WHERE SPECIFIED.
- f. 1/2" DIA. FENDER WASHER SHALL HAVE A MINIMUM OUTSIDE DIAMETER OF 1 1/2" AND MINIMUM THICKNESS OF 7/64".
- g. ALL EXPOSED STEEL FINISH SHALL BE DETERMINED BY THE OWNER/END USER. IF UNKNOWN, PROVIDE STAINLESS STEEL.

6. STEEL STUDS, JOISTS AND ACCESSORIES:

- a. STEEL STUDS, JOISTS AND ACCESSORIES SHALL BE MANUFACTURED BY MEMBERS OF STEEL STUD MANUFACTURER'S ASSOCIATION (SSMA ICC-ES ESR-3064P PRODUCT TECHNICAL GUIDE).
- b. STEEL STUDS, JOISTS AND ACCESSORIES SHALL BE FORMED FROM STEEL WITH A MINIMUM YIELD STRESS OF 50 KSI FOR 16 GAUGE AND HEAVIER ITEMS, OR 33 KSI FOR 18 GAUGE AND LIGHTER ITEMS.
- c. STEEL STUDS, JOISTS AND ACCESSORIES SHALL BE GALVANIZED IN ACCORDANCE WITH ASTM A653, G60, UNLESS OTHERWISE NOTED.
- d. ALL WELDS OF GALVANIZED STEEL STUDS, JOISTS AND ACCESSORIES SHALL BE TOUCHED UP WITH A ZINC RICH PAINT.
- e. SHEET METAL SCREWS (SMS) SHALL COMPLY WITH THE REQUIREMENTS OF STEEL STUD MANUFACTURER'S ASSOCIATION (SSMA) ICC-ES ESR-3064P AND MAY BE SELF-DRILLING AND/OR SELF-TAPPING AS DESIRED.
- f. PENETRATION OF SCREWS THROUGH JOINED MATERIALS SHALL NOT BE LESS THAN 3 EXPOSED THREADS. SCREWS SHALL BE INSTALLED AND TIGHTENED IN ACCORDANCE WITH SCREW MANUFACTURER'S RECOMMENDATIONS.
- g. SHEET METAL WELDING SHALL BE PERFORMED IN ACCORDANCE WITH AWS D1.3. ELECTRODES SHALL CONFORM TO AWS A5.1 OR A5.5 AND TABLE 5.2 OF AWS D1.3.



SHEET TITLE: GENERAL NOTES		
EQUIPMENT INFORMATION: MODELS EC770, EC770L, EC1260BL, EC1540, EC1540BL, EC1730BL, EC2060, EC2140 BL		
 1 Peters Canyon Rd. Irvine, CA 92606 www.idsgi.com 949.387.8500		ISSUE DATE: JUNE 6, 2014
		OPM SHEET NO: 1.2

ENG STAMP:


7. WOOD:

- a. WOOD MEMBERS SHALL BE DOUGLAS FIR-LARCH PER WCLIB VISUALLY GRADED DIMENSION LUMBER AND SHALL BE SURFACED DRY (19% MOISTURE CONTENT MAXIMUM). STRUCTURAL FRAMING MEMBERS SHALL BE S4S AND GRADE MARKED AS No. 2 MINIMUM.
- b. NAILS SHALL BE COMMON WIRE NAILS (0.131" DIA. x 2 1/2" FOR 8d; 0.148" DIA. x 3" FOR 10d; 0.148" DIA. x 3 1/4" FOR 12d; 0.162" DIA. x 3 1/2" FOR 16d) OR FASTENERS PROVIDED WITH HARDWARE CONNECTORS SPECIFIED IN OPM.
- c. A34 CONNECTOR AND ACCESSORIES SHALL BE SIMPSON STRONG-TIE (IAPMO ER-112).
- d. WOOD SCREWS SHALL BE SCREWED, NOT DRIVEN, INTO WOOD MEMBER. WOOD SCREW SHALL BE SIMPSON STRONG-TIE SDS SCREW 1/4" DIA. x 3 1/2" MIN (ICC-ES ESR-2236) WITH 2 1/4" MIN PENETRATION INTO WOOD MEMBER.
- e. WOOD SCREW FINISH SHALL BE DETERMINED BY THE OWNER/END USER. IF UNKNOWN, PROVIDE STAINLESS STEEL.

8. POST-INSTALLED EXPANSION ANCHORS:

- a. FOR ANCHORAGE TO CMU: ALLOWABLE STRENGTHS SHALL BE COMPARED TO ALLOWABLE STRESS DESIGN (ASD) LEVEL DEMAND IN ACCORDANCE WITH ASCE 7-10 SECTION 12.4 AND 2013 CALIFORNIA BUILDING CODE SECTIONS 1908A.

FOR ANCHORAGE TO CONCRETE: STRENGTH DESIGN IN ACCORDANCE WITH ASCE 7-10 SECTION 12.4 AND 2013 CALIFORNIA BUILDING CODE SECTIONS 1909A.
- b. SHALL HAVE SPECIAL INSPECTION AND TESTING IN ACCORDANCE WITH THE 2013 CALIFORNIA BUILDING CODE SECTIONS 1705A.3 & 1913A.7.3.
- c. MINIMUM CONCRETE STRENGTH $f_c = 3000$ PSI OR MINIMUM SOLID GROUTED CMU STRENGTH $f_m = 1500$ PSI. SEOR SHALL PROVIDE CALCULATIONS JUSTIFYING THE CMU WALL REMAINS ELASTIC.
- d. POST-INSTALLED ANCHORS SHALL NOT BE USED IN PRE-STRESSED CONCRETE UNLESS NON-DESTRUCTIVE TESTING METHODS ARE USED TO LOCATE STRAND AND REINFORCING PRIOR TO ANCHOR INSTALLATION.
- e. POST-INSTALLED ANCHOR INSTALLATION SHALL NOT NICK OR DAMAGE EXISTING REINFORCEMENT. SHOULD THIS OCCUR THE SEOR BE NOTIFIED IMMEDIATELY. EXPANSION ANCHORS SHALL BE INSTALLED 1" CLEAR OF EXISTING REINFORCEMENT.
- f. POST-INSTALLED ANCHORS SHALL BE INSTALLED PER TABLES 1.5A AND 1.5B ON SHEET 1.5.
- g. ALL POST-INSTALLED ANCHOR VALUES ARE FOR CRACKED CONCRETE.
- h. POST-INSTALLED ANCHOR FINISH SHALL BE DETERMINED BY THE END USER. IF UNKNOWN, PROVIDE STAINLESS STEEL.



<p>ENG STAMP:</p>	SHEET TITLE: GENERAL NOTES	
	EQUIPMENT INFORMATION: MODELS EC770, EC770L, EC1260BL, EC1540, EC1540BL, EC1730BL, EC2060, EC2140 BL	
	<p>1 Peters Canyon Rd. Irvine, CA 92606 www.idsgi.com 949.387.8500</p>	
		ISSUE DATE: JUNE 6, 2014
		OPM SHEET NO: 1.3

POST INSTALLED EXPANSION ANCHORS CONT.:

- i. CONCRETE ANCHORS AND ICC-ES REPORTS:
 HILTI KWIK BOLT TZ (KBTZ): ESR-1917
 SIMPSON STRONG-TIE STRONG BOLT 2 (SB2): ESR-3037
 POWERS POWER-STUD+ SD1: ESR-2818
 POWERS POWER-STUD+ SD2: ESR-2502

- j. CONCRETE MASONRY UNIT (CMU) ANCHORS AND ICC-ES REPORTS:
 HILTI KWIK BOLT 3 (KB3): ESR-1385



- k. ANCHOR EMBEDMENT DEPTHS SHOWN IN TABLE 1.6b THROUGH 1.6b ARE THE MINIMUM REQUIRED TO MEET THE SPECIFIC LOADING REQUIREMENTS. SEE DIAGRAM BELOW FOR NOMENCLATURE.

- l. POST-INSTALLED ANCHORS SHALL BE TESTED PER 2013 CBC SECTION 1913A.7 TO TORQUE OR DIRECT-PULL TENSION VALUES NOTED FOR THE ANCHOR INSTALLED.

- m. ALL (100%) OF THE POST-INSTALLED ANCHORS SHALL BE TESTED. TESTING SHALL BE WITNESSED AND REPORTED BY THE INSPECTOR OF RECORD OR SPECIAL INSPECTOR. REPORTS OF TEST RESULTS TO BE SUBMITTED TO OSHPD.

- n. TEST ACCEPTANCE IS BY REACHING THE TORQUE WITHIN "1/2 TURN OF THE NUT" OR MAINTAINING THE TEST LOAD FOR AT LEAST 15 SECONDS AND SHALL EXHIBIT NO DISCERNABLE MOVEMENT.

- o. IF STAINLESS STEEL ANCHOR IS REQUIRED BY THE OWNER/END USER, UTILIZE MANUFACTURER'S SAME PRODUCT WITH TENSION TEST UNLESS INCLUDED IN ICC-ES REPORT.



SHEET TITLE:
GENERAL NOTES

EQUIPMENT INFORMATION:
MODELS EC770, EC770L, EC1260BL, EC1540, EC1540BL, EC1730BL, EC2060, EC2140 BL

IDS GROUP
1 Peters Canyon Rd.
Irvine, CA 92606
www.idsgi.com 949.387.8500

ENTHERMICS
Medical Systems

ISSUE DATE:
JUNE 6, 2014

OPM SHEET NO:
1.4

ENG STAMP:

TABLE 1.5a

MINIMUM REQUIRED ANCHOR OPTIONS FOR CONCRETE WALL (f'c = 3,000 PSI MIN AND hmin PER TABLE BELOW) SEE SHEET 1.4 NOTE 8i FOR RELATED RESEARCH REPORT								
ANCHOR SIZE	hef	hnom	ho	hmin	min edge distance	min spacing	INSTALLATION AND TESTING TORQUE	TENSION TEST LOAD REQUIREMENT
1/2" DIA KBTZ	3 1/4"	3 5/8"	4"	6"	4"	8"	40 ft-lb	3750 LB AT NORMAL WEIGHT CONCRETE
1/2" DIA SB2	3 3/8"	3 7/8"	4 1/8"	6"	4"	8"	60 ft-lb	
1/2" DIA SD2	3 1/4"	3 3/4"	4"	5 3/4"	4"	8"	40 ft-lb	
1/2" DIA SD1	3 1/4"	3 3/4"	4"	6"	4"	8"	40 ft-lb	

TABLE 1.5b

MINIMUM REQUIRED ANCHOR OPTIONS FOR SOLID GROUTED CMU WALL (f'm = 1,500 PSI MIN AND hmin PER TABLE BELOW) SEE SHEET 1.4 NOTE 8j FOR RELATED RESEARCH REPORT							
ANCHOR SIZE	hnom	ho	hmin	min edge distance	min spacing	INSTALLATION AND TESTING TORQUE	TENSION TEST LOAD REQUIREMENT
1/2" DIA KB3	3 1/2"	4"	6" nom.	4"	8"	25 ft-lb	1100 LB

NOTES:

1. EXPANSION ANCHORS DESIGNED TO ICC-ES AC01 ARE LIMITED TO ALLOWABLE STRESS DESIGN **ONLY**, IN ACCORDANCE WITH AC01 1.2. HENCE, STRENGTH DESIGN VALUES ARE NOT ACCEPTABLE. ALLOWABLE STRESS VALUES CAN BE SHOWN PROVIDED, SEOR WILL VERIFY THAT:
 - A) MASONRY IS NOT CRACKED AS DEFINED IN ICC-ES AC01 SECTION 2.3. SEOR SHALL PROVIDE CALCULATIONS TO SHOW THAT THE MASONRY WALL WOULD NOT CRACK UNDER THE DESIGN EARTHQUAKE LOADS UNDER ALL SERVICE LOAD CONDITIONS; THE WALL HAS TO REMAIN ELASTIC;
 - B) MASONRY IS FULLY GROUTED IN ACCORDANCE WITH ESR-1385 SECTION 3.2;
 - C) CONDITION OF USE REQUIREMENTS IN ACCORDANCE WITH ESR-1385 SECTION 5.0 IS SATISFIED.



ENG STAMP:

SHEET TITLE:
GENERAL NOTES

EQUIPMENT INFORMATION:
MODELS EC770, EC770L, EC1260BL, EC1540, EC1540BL, EC1730BL, EC2060, EC2140 BL



ISSUE DATE:
JUNE 6, 2014

OPM SHEET NO:
1.5


MODEL	Cabinet Weight (pounds)	Blanket Space (Cubic feet per MRF)	Blanket Weight (pounds)	Fluid Capacity (liters)	Fluid Weight (pounds)	Content Weight (pounds)	Wp (pounds)
EC770	178	7.7	154	n/a	n/a	154	332
EC770L	250	n/a	n/a	72	159	159	409
EC1260BL	433	8.7	174	18	40	214	647
EC1540	412	15.4	308	n/a	n/a	308	720
EC1540BL	414	7.7	154	72	159	313	727
EC1730BL	360	11.9	238	38	84	322	682
EC2060	356	20.7	414	n/a	n/a	414	770
EC2140BL	441	15.5	310	25	55	365	806

NOTES:

- SEE SHEET 2.3 FOR NUMBER OF ATTACHMENT BRACKETS.
- SEE SHEET 2.4 AND 2.5 FOR ADJACENT SUPPORT STRUCTURE REFERENCE FOR ATTACHMENT DETAILING.
- SEE SHEETS 2.6 AND 2.7 SHOWING REFERENCE MATERIAL FOR CUSTOM FABRICATED ATTACHMENT BRACKET ASSEMBLIES SUPPLIED BY ENTHERMICS:
 TYPE A: ENTHERMICS MODEL 5017207 (OSHPD ASSEMBLY 2.5")
 TYPE B: ENTHERMICS MODEL 5017226 (OSHPD ASSEMBLY 5.0")
- EACH CUSTOM FABRICATED ATTACHMENT BRACKET TYPE A OR B IS ATTACHED TO THE 20 GAGE MIN CABINET FRAMING/SHELL WITH ELEVEN (11) # 10 SHEET METAL SCREWS (9 SCREWS TO THE BACK OF THE CABINET, AND 2 SCREWS TO THE SIDE OF THE CABINET) SEE SHEETS 8.1 AND 8.2 FOR ADDITIONAL INFORMATION.
- CUSTOM FABRICATED ATTACHMENT BRACKET TYPE B HAS SIX (6) #14 SMS FROM THE BOTTOM OF PART (b) TO THE SUPPORT TABS OF PART (d). ALSO, SIX (6) #10 SMS FROM THE SIDE TABS OF PART (d) TO THE SIDES OF PART (b).



SHEET TITLE: CABINET CONTENTS AND WEIGHTS	
EQUIPMENT INFORMATION: MODELS EC770, EC770L, EC1260BL, EC1540, EC1540BL, EC1730BL, EC2060, EC2140 BL	
 1 Peters Canyon Rd. Irvine, CA 92606 www.idsgi.com 949.387.8500	 ISSUE DATE: JUNE 6, 2014 OPM SHEET NO: 2.1



TYPICAL CABINET ISOMETRIC

MODEL	Width (inches)	Depth (inches)	Height (inches)	"B" (inches)	"T" (inches)	"S" (inches)	BRACKET	TOP BRACKET C.G., Z (in)	BOTTOM BRACKET C.G., Z (in)
EC770	22.2	28.8	40.3	12 ± 3	6 ± 3	2 1/2 ±	TYPE A	26.6	20.2
EC770L	22.2	28.0	40.7	12 ± 3	6 ± 3	5 ±	TYPE B	26.9	20.6
EC1260BL	30.1	20.2	69.2	12 ± 3	12 ± 3	2 1/2 ±	TYPE A	45.7	34.6
EC1540	23.9	28.8	72.7	12 ± 3	6 ± 3	2 1/2 ±	TYPE A	48.0	36.4
EC1540BL	22.1	28.8	72.5	12 ± 3	6 ± 3	5 ±	TYPE B	47.8	36.2
EC1730BL	30.2	25.4	70.9	14 ± 3	8 ± 3	2 1/2 ±	TYPE A	46.8	35.4
EC2060	30.2	25.2	70.8	14 ± 3	8 ± 3	2 1/2 ±	TYPE A	47.2	35.4
EC2140BL	30.0	26.0	73.5	14 ± 3	8 ± 3	2 1/2 ±	TYPE A	48.9	36.7



SHEET TITLE:

CABINET CONFIGURATION AND DIMENSIONS

EQUIPMENT INFORMATION:

MODELS EC770, EC770L, EC1260BL, EC1540, EC1540BL, EC1730BL, EC2060, EC2140 BL



ISSUE DATE:

JUNE 6, 2014

OPM SHEET NO.:

2.2

ENG STAMP:

NOTES:

1. MODELS EC770 AND EC770L SHALL BE ATTACHED WITH 4 BRACKETS PER SHEET 2.4 FOR ALL Z/H \leq 1.0 AND ANY Sds \leq 2.5

2. MODELS EC1260BL, EC1540, EC1540BL, EC1730BL, EC2060, AND EC2140BL SHALL BE ATTACHED WITH 4 BRACKETS (SHEET 2.4) EITHER 6 OR 8 BRACKETS (SHEET 2.5) PER THE TABLE BELOW.

NUMBER OF ATTACHMENT BRACKETS						
Z/H RATIO	Sds VALUES LESS THAN OR EQUAL TO 2.5					
	<1.4	\leq 1.5	\leq 1.75	\leq 2.0	\leq 2.25	\leq 2.5
\leq 0.4	4	4	4	4	4	6
\leq 0.5	4	4	4	4	6	6
\leq 0.6	4	4	4	6	6	SEE NOTE 4
\leq 0.8	4	4	6	6		
0.8 < x \leq 1.0	4	6	6			

3. Z = HEIGHT IN STRUCTURE OF POINT OF ATTACHMENT OF COMPONENT WITH RESPECT TO THE BASE. FOR ITEMS AT OR BELOW THE BASE, Z SHALL BE TAKEN AS 0. THE VALUE OF Z/H NEED NOT EXCEED 1.0


H = AVERAGE ROOF HEIGHT OF STRUCTURE WITH RESPECT TO THE BASE.

4. PROVIDE 6 CONNECTIONS (2 AT THE TOP AND 1 AT THE BOTTOM) TO NON CONCRETE OR MASONRY WALLS.


PROVIDE 8 CONNECTIONS (2 AT THE TOP AND 2 AT THE BOTTOM) TO CONCRETE OR MASONRY WALLS.



SHEET TITLE: NUMBER OF ATTACHMENT BRACKETS REQUIRED		
EQUIPMENT INFORMATION: MODELS EC770, EC770L, EC1260BL, EC1540, EC1540BL, EC1730BL, EC2060, EC2140 BL		
<p>IDS GROUP 1 Peters Canyon Rd. Irvine, CA 92606 www.idsgi.com 949.387.8500</p>	<p>ENTHERMICS Medical Systems</p>	ISSUE DATE: JUNE 6, 2014
		OPM SHEET NO: 2.3



TYPICAL CABINET PLAN VIEW



TYPICAL CABINET SIDE ELEVATION

NOTE:

- "Y" = 2.6" MIN, 6.9" MAX
- MAXIMUM FORCES SHOWN ARE TO THE FACE OF ARCHITECTURAL FINISH. THE THICKNESS OF THE ARCHITECTURAL FINISH SHOULD BE TAKEN INTO ACCOUNT WHEN DESIGNING WALL STUDS.



SHEET TITLE:
4 ATTACHMENTS - CABINET TOP VIEW AND SIDE ELEVATION


EQUIPMENT INFORMATION:
 MODELS EC770, EC770L, EC1260BL, EC1540, EC1540BL, EC1730BL, EC2060, EC2140 BL

ISSUE DATE:
 JUNE 6, 2014


OPM SHEET NO:
2.4

IDS GROUP
 1 Peters Canyon Rd.
 Irvine, CA 92606
 www.idsgi.com 949.387.8500

ENTHERMICS
 Medical Systems



TYPICAL CABINET PLAN VIEW



TYPICAL CABINET SIDE ELEVATION

NOTE:

- "Y" = 2.6" MIN, 6.9" MAX
- MAXIMUM FORCES SHOWN ARE TO THE FACE OF ARCHITECTURAL FINISH. THE THICKNESS OF THE ARCHITECTURAL FINISH SHOULD BE TAKEN INTO ACCOUNT WHEN DESIGNING WALL STUDS.



ENG STAMP:

SHEET TITLE:

6 OR 8 ATTACHMENTS - CABINET TOP VIEW AND SIDE VIEW

EQUIPMENT INFORMATION:

MODELS EC770, EC770L, EC1260BL, EC1540, EC1540BL, EC1730BL, EC2060, EC2140 BL



ISSUE DATE:

JUNE 6, 2014

OPM SHEET NO.:

2.5


NOTES:

1. GENERAL ASSEMBLY INFORMATION PROVIDED. SPECIFIC PRODUCT INFORMATION AND INSTALLATION INSTRUCTIONS PER ENTHERMICS MEDICAL SYSTEMS.

2. SEE SHEETS 3.2, 4.2, OR 7.2 FOR HOLES UTILIZED FOR ATTACHMENT.

3. SEE SHEET 2.5 NOTES 4 AND 5 FOR REQUIRED CONNECTIONS BETWEEN ASSEMBLY.

- (a) STAINLESS STEEL 16GA BRACKET WITH 11GA GUSSET (SPOT WELDED). FIELD ATTACHED TO CABINET WITH PROVIDED SMS.
- (b) STAINLESS STEEL 16GA BRACKET WITH 11GA GUSSET (SPOT WELDED) FIELD ATTACHMENT TO SUPPORT STRUCTURE WITH INSTALLER PROVIDED HARDWARE.
- (c) 1/2" DIAMETER STAINLESS STEEL QUICK RELEASE PIN, SEE PAGE 1.2 NOTE 5a.




ASSEMBLY ISOMETRIC

PROVIDE 3/4" LONG 16 GA FILLET WELD EACH SIDE OF BRACKET TO STIFFENER PLATE, TYP OR SCREW CONNECTION SHOWN.

PROVIDE (2) 8-32 MACHINE SCREWS THROUGH BRACKET AND INTO STIFFENERS EACH SIDE, TYP OR FILLET WELDS SHOWN. DRILL AND TAP STIFFENER PLATE.

SEE SHEET 8.1 FOR ATTACHMENT

SEE SHEET 8.2 FOR ATTACHMENT



EXPLODED VIEW



ENG STAMP:

SHEET TITLE:
REFERENCE BRACKET ASSEMBLY TYPE A

EQUIPMENT INFORMATION:
MODELS EC770, EC770L, EC1260BL, EC1540, EC1540BL, EC1730BL, EC2060, EC2140 BL

IDS GROUP
1 Peters Canyon Rd.
Irvine, CA 92606
www.idsgi.com 949.387.8500

ENTHERMICS
Medical Systems

ISSUE DATE:
JUNE 6, 2014


OPM SHEET NO:
2.6

NOTES:

1. GENERAL ASSEMBLY INFORMATION PROVIDED. SPECIFIC PRODUCT INFORMATION AND INSTALLATION INSTRUCTIONS PER ENTHERMICS MEDICAL SYSTEMS.

2. SEE SHEETS 3.2, 4.2, OR 7.2 FOR HOLES UTILIZED FOR ATTACHMENT.


3. SEE SHEET 2.5 NOTES 4 AND 5 FOR REQUIRED CONNECTIONS BETWEEN ASSEMBLY.



ASSEMBLY ISOMETRIC

PROVIDE (2) 8-32 MACHINE SCREWS THROUGH BRACKET AND INTO STIFFENERS EACH SIDE, TYP OR FILLET WELDS SHOWN. DRILL AND TAP STIFFENER PLATE.

- (a) STAINLESS STEEL 16 GA BRACKET WITH 11 GA GUSSET (SPOT WELDED). FIELD ATTACHED TO CABINET WITH PROVIDED SMS.
- (b) STAINLESS STEEL 16 GA BRACKET WITH 11 GA GUSSET (SPOT WELDED) FIELD ATTACHED TO PART D WITH PROVIDED SMS.
- (c) 1/2" DIAMETER STAINLESS STEEL QUICK RELEASE PIN, SEE PAGE 1.2 NOTE 5a.
- (d) STAINLESS STEEL 16 GA EXTENSION. FIELD ATTACHED TO SUPPORT STRUCTURE WITH INSTALLER PROVIDED HARDWARE.



EXPLODED VIEW

SEE SHEET 8.1 FOR ATTACHMENT

SEE SHEET 8.2 FOR ATTACHMENT




SHEET TITLE: **REFERENCE BRACKET ASSEMBLY TYPE B**

EQUIPMENT INFORMATION: **MODELS EC770, EC770L, EC1260BL, EC1540, EC1540BL, EC1730BL, EC2060, EC2140 BL**




ISSUE DATE: **JUNE 6, 2014**

OPM SHEET NO: **2.7**



16 GA. MIN. METAL STUD WALL WITH 1 OR 2 LAYERS OF 5/8" ARCHITECTURAL FINISH BY OTHERS

BRACKET TYPE A



BACKING PLATE/STUD SEE SHEET 3.3

BRACKET TYPE B



ENG STAMP:


SHEET TITLE:
ATTACHMENT TO METAL STUD WALL

EQUIPMENT INFORMATION:
MODELS EC770, EC770L, EC1260BL, EC1540, EC1540BL, EC1730BL, EC2060, EC2140 BL




ISSUE DATE:
JUNE 6, 2014

OPM SHEET NO:
3.1



BRACKET TYPE A



BRACKET TYPE B

KEY NOTES:

- 1** # 14 (1/4" DIA.) SHEET METAL SCREW (SMS)
- 2** 9/16" DIA. HOLE NOT USED



ENG STAMP:

SHEET TITLE:
BRACKET TO SHEET METAL


EQUIPMENT INFORMATION:
MODELS EC770, EC770L, EC1260BL, EC1540, EC1540BL, EC1730BL, EC2060, EC2140 BL

IDS GROUP
1 Peters Canyon Rd.
Irvine, CA 92606
www.idsgj.com 949.387.8500


ENTHERMICS
Medical Systems

ISSUE DATE:
JUNE 6, 2014

OPM SHEET NO:
3.2



ISOMETRIC VIEW



PLAN VIEW

NOTES

- 1. NOTCHING OR CUTTING OF BACKING IS NOT PERMITTED
- 2. DO NOT ATTACH TWO BRACKETS IN SPACE BETWEEN STUDS (ONE STUD SPACING IS MINIMUM SPACING BETWEEN BRACKETS ON ADJACENT UNITS).




SHEET TITLE: **BACKING TRACK CONFIGURATION**

EQUIPMENT INFORMATION: **MODELS EC770, EC770L, EC1260BL, EC1540, EC1540BL, EC1730BL, EC2060, EC2140 BL**




ISSUE DATE: **JUNE 6, 2014**

OPM SHEET NO: **3.3**



BRACKET TYPE A

$f_c = 3000$ PSI MIN CONCRETE
SUPPORT WALL BY OTHERS



BRACKET TYPE B



ENG STAMP:

SHEET TITLE:

ATTACHMENT TO CONCRETE WALL

EQUIPMENT INFORMATION:

MODELS EC770, EC770L, EC1260BL, EC1540, EC1540BL, EC1730BL, EC2060, EC2140 BL




ISSUE DATE:

JUNE 6, 2014

OPM SHEET NO:


4.1



KEY NOTES:

- 1** 1/2" DIA. EXP. BOLT WITH FENDER WASHER (SEE SHEET 1.5 TABLES 1.5a AND 1.5b AS APPLICABLE)
- 2** 1/4" DIA. CONC. SCREW (NON-STRUC)
- 3** 0.273" MAX DIA. HOLE NOT USED

BRACKET TYPE A



BRACKET TYPE B


NOTES:

1. SEE SHEET 3.2 FOR BRACKET HOLE LOCATION DIMENSIONS AND MATERIAL TYPE.




ENG STAMP:

SHEET TITLE: BRACKET TO CEMENTITIOUS MATERIAL	
EQUIPMENT INFORMATION: MODELS EC770, EC770L, EC1260BL, EC1540, EC1540BL, EC1730BL, EC2060, EC2140 BL	
<p>IDS GROUP 1 Peters Canyon Rd. Irvine, CA 92606 www.idsgi.com 949.387.8500</p>	<p>ENTHERMICS Medical Systems</p>
	<p>ISSUE DATE: JUNE 6, 2014</p> <p>OPM SHEET NO: 4.2</p>



BRACKET TYPE A

f'm = 1500 PSI MIN SOLID
GROUTED CONCRETE
MASONRY UNIT (CMU)
SUPPORT WALL BY OTHERS



BRACKET TYPE B



SHEET TITLE:
ATTACHMENT TO CMU WALL


EQUIPMENT INFORMATION:
MODELS EC770, EC770L, EC1260BL, EC1540, EC1540BL, EC1730BL, EC2060, EC2140 BL




ISSUE DATE:
JUNE 6, 2014

OPM SHEET NO:
5.1

ENG STAMP:



BRACKET TYPE A



BRACKET TYPE B

NOTE:
 1. SCREWS SHOWN ARE SCREWS TO THE BACKING PLATE. THE SCREWS ATTACHING THE BACKING PLATE TO THE METAL STUDS ARE NOT SHOWN FOR CLARITY. SEE SHEET 6.2 FOR SCREW REQUIRED TO ATTACH BACKING PLATE TO WALL STUDS.




SHEET TITLE:
BRACKET TO SURFACE MOUNTED PLATE

EQUIPMENT INFORMATION:
 MODELS EC770, EC770L, EC1260BL, EC1540, EC1540BL, EC1730BL, EC2060, EC2140 BL




ISSUE DATE:
 JUNE 6, 2014

OPM SHEET NO:
6.1



NOTES

1. NOTCHING OR CUTTING OF BACKING PLATE IS NOT PERMITTED
2. BACKING PLATE IS PLACED ON THE SURFACE OF 1 OR 2 LAYERS OF 5/8" ARCHITECTURAL FINISH.
3. DO NOT ATTACH TWO BRACKETS IN SPACE BETWEEN STUDS (ONE STUD SPACING IS THE MINIMUM SPACING BETWEEN BRCKETS ON ADJACENT UNITS).



SHEET TITLE: SURFACE PLATE CONFIGURATION		
EQUIPMENT INFORMATION: MODELS EC770, EC770L, EC1260BL, EC1540, EC1540BL, EC1730BL, EC2060, EC2140 BL		
 IDS GROUP 1 Peters Canyon Rd. Irvine, CA 92606 www.idsgi.com 949.387.8500	 ENTHERMICS Medical Systems	ISSUE DATE: JUNE 6, 2014
		OPM SHEET NO: 6.2

SIMPSON A34, TYP
TOP AND BOTTOM
EACH END

BACKING SEE
SHEET 7.3

2 1/2" +/-

7.2

8.1

ASSEMBLY BRACKET

CABINET (NOT
IN OPM SCOPE)

8.2

WOOD STUD WALL WITH
1 OR 2 LAYERS OF 5/8"
ARCHITECTURAL FINISH
BY OTHERS

BRACKET TYPE A

SIMPSON A34, TYP
TOP AND BOTTOM
EACH END

BACKING SEE
SHEET 7.3

5" +/-

7.2

8.1

ASSEMBLY BRACKET

CABINET (NOT
IN OPM SCOPE)

8.2

BRACKET TYPE B

SHEET TITLE:

ATTACHMENT TO WOOD STUD WALL

EQUIPMENT INFORMATION:

MODELS EC770, EC770L, EC1260BL, EC1540, EC1540BL, EC1730BL, EC2060, EC2140 BL

ISSUE DATE:

JUNE 6, 2014

OPM SHEET NO:

7.1




ENG STAMP:




KEY NOTES:

- 1** 1/4" DIA. WOOD SCREW (0.273" MAX HOLE IN STEEL MEMBER PER ICC-ES ESR REPORT)
- 2** 9/16" DIA. HOLE NOT USED



BRACKET TYPE A




BRACKET TYPE B

TAB ABOVE, TYP.

NOTES:


1. SEE SHEET 3.2 FOR BRACKET HOLE LOCATION DIMENSIONS & MATERIAL TYPE.




ENG STAMP:

SHEET TITLE: BRACKET TO WOOD	
EQUIPMENT INFORMATION: MODELS EC770, EC770L, EC1260BL, EC1540, EC1540BL, EC1730BL, EC2060, EC2140 BL	
 <p>IDS GROUP 1 Peters Canyon Rd. Irvine, CA 92606 www.idsgi.com 949.387.8500</p>	 <p>ENTHERMICS Medical Systems</p>
	ISSUE DATE: JUNE 6, 2014 OPM SHEET NO: 7.2


2x CRIPPLE STUD REQ. AT (E) 2X AND 3X STUD FRAMING WHERE A34 CLIPS ARE PROVIDED ON EACH SIDE OF STUD.



ISOMETRIC DETAIL



CONTINUOUS BACKING BETWEEN (3) STUDS MIN AND TO THE NEXT STUD PAST EA. END OF THE RESTRAINED CABINET.




PLAN VIEW

NOTES:

- 1. DO NOT ATTACH TWO BRACKETS IN SPACE BETWEEN STUDS (ONE STUD SPACING IS MINIMUM SPACING BETWEEN BRACKETS ON ADJACENT UNITS).



SHEET TITLE: WOOD BACKING CONFIGURATION	
EQUIPMENT INFORMATION: MODELS EC770, EC770L, EC1260BL, EC1540, EC1540BL, EC1730BL, EC2060, EC2140 BL	
<p>IDS GROUP 1 Peters Canyon Rd. Irvine, CA 92606 www.idsgj.com 949.387.8500</p>	<p>ENTHERMICS Medical Systems</p>
	<p>ISSUE DATE: JUNE 6, 2014</p> <p>OPM SHEET NO: 7.3</p>



ENG STAMP:

SHEET TITLE:

BRACKET TYPE A OR B ATTACHMENT TO BACK OF CABINET

EQUIPMENT INFORMATION:

MODELS EC770, EC770L, EC1260BL, EC1540, EC1540BL, EC1730BL, EC2060, EC2140 BL




ISSUE DATE:

JUNE 6, 2014

OPM SHEET NO:

8.1



NOTE:
BRACKET PART (b) IS NOT SHOWN FOR CLARITY.



SHEET TITLE:
BRACKET TYPE A OR B ATTACHMENT TO SIDE OF CABINET

EQUIPMENT INFORMATION:
MODELS EC770, EC770L, EC1260BL, EC1540, EC1540BL, EC1730BL, EC2060, EC2140 BL



ISSUE DATE:
JUNE 6, 2014

OPM SHEET NO:
8.2