

OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT FACILITIES DEVELOPMENT DIVISION

OFFICE USE ONLY APPLICATION FOR OSHPD PREAPPROVAL OF **MANUFACTURER'S CERTIFICATION (OPM) APPLICATION #:** OPM-0345-13

OSHPD Preapproval of Manufacturer's Certification (OPM)
Type: ☐ New ☐ Renewal ☐ Update to Pre-CBC 2013 OPA Number: 2498-07
Manufacturer Information
Manufacturer: BD Life Sciences – Diagnostic Systems
Manufacturer's Technical Representative: David Lentz
Mailing Address: 7 Loveton Circle – mc 622, Sparks, MD 21152
Telephone: 410-316-4779 Email: Dave_lentz@bd.com
Product Information
Product Name: BD BACTEC™ MGIT™ 320
Product Type: Mycobacteria Culture System Unit mounted on stand
Product Model Number: BD BACTEC™ MGIT™ 320
General Description: Tuberculosis Detection Instrument
DATE: 07/14/2016
Applicant Information
Applicant Company Name: BD Life Sciences - Diagnostic Systems
Contact Person: David Lentz
Mailing Address: 7 Loveton Circle – mc 622, Sparks, MD 21152
Telephone: 410-316-4779 Email: Dave_lentz@bd.com
I hereby agree to reimburse the Office of Statewide Health Planning and Development review fees in accordance with the California Administrative Code, 2013.
Signature of Applicant: Date: 05/02/2016
Title: Director, R&D Hardware Engineering Company Name: BD Life Sciences – Diagnostic Systems

"Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"







OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT FACILITIES DEVELOPMENT DIVISION

Registered Design Professional Preparing Engineering Recommendations								
Company Name: _C	CYS Structural Engineers, Inc.							
Name: D	Dieter T. Siebald	Califorr	nia License Number:	S4346				
Mailing Address: 2	2495 Natomas Park Drive, Suite #650,	Sacramento	o, CA 95833					
Telephone:	916-920-2020	Email: d	dieters@cyseng.com					
OSHPD Special Se	eismic Certification Preapproval ((OSP)						
(Separate applic	Certification is preapproved under OSI cation for OSP is required) Certification is not preapproved	P-						
Certification Metho	od(s)	CODE						
☐ Testing in accord	dance with:	FM 195	50-10					
supports and attachm bracings, test criteria of Analysis Experience Data	than those adopted by the California Branch are not permitted. For distribution other than those adopted in the CBSC DATE: 07/ Testing, Analysis, and/or Experience Date.	system, into 2013 may b	erior partitio <mark>n w</mark> all, and s be used when approved	suspended ceiling seismic				
List of Attachment	ts Supporting the Manufacturer's	Certificat	ion					
☐ Test Report ☐ Other(s) (Plea	☑ Drawings ☑ Calculase Specify):	lations	☐ Manufacturer's Ca	atalog				
OFFICE USE ONLY -	OSHPD APPROVAL VALID FOR C	BC 2013 ON	ILY					
Signature:	ey Kikumoto		Date:	07-14-2016				
supports and attachm bracings, test criteria of the composition of Total Combination of Total	Testing, Analysis, and/or Experience Dates Supporting the Manufacturer's Drawings Calculates Specify): OSHPD APPROVAL VALID FOR CE	ata (Please Certifications	erior partition wall, and soe used when approved ase Specify): Manufacturer's Ca	suspended ceiling seismic by OSHPD prior to testing.				

"Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"





Page 2 of 2



TABLE OF CONTENTS OPM-0345-13

	PAGE
GENERAL NOTES	2
SYSTEM OVERVIEW & DESIGN CRITERIA	. 4
ABBREVIATIONS	5
STAND CONFIGURATION ELEVATIONS	6
ATTACHMENT PLANS UNIT TO STAND STAND TO FLOOR	7 8
BRACKET HEAVY ANCHOR DETAIL FRONT MOUNT HEAVY ANCHOR DETAIL STAND ASSEMBLY DETAIL LEVELING LEG DETAILS	9 9 10 11
ATTACHMENT DETAILS STAND TO CONCRETE FILL OVER METAL DECK (CASE 1) STAND TO 4" SLAB ON GRADE (CASE 2)	. 12 14

NOTES: 1. THESE DRAWINGS ARE PREPARED FOR BUILIFE SCIENCES - DIAGNOSTIC SYSTEMS, SPARKS, MARYLAND.

- 2. THE CONTRACTOR AND INSPECTOR OF RECORD SHALL OBTAIN A COPY OF THIS PRE-APPROVAL FROM THE OFFICE OF STATEWIDE HEALTH PLANNING & DEVELOPMENT (OSHPD) PRE-APPROVAL PROGRAMS WEBSITE.
- 3. THIS PRE-APPROVAL COVERS THE SUPPORTS AND ATTACHMENTS OF THE EQUIPMENT TO THE SUPPORTING STRUCTURE. THE EQUIPMENT, STAND & ATTACHMENT HARDWARE ARE SUPPLIED BY THE MANUFACTURER. THE EXPANSION ANCHORS, THRU-BOLTS & STRUT PLATES SHOWN IN THIS OPM SHALL BE SUPPLIED AND INSTALLED BY THE CONTRACTOR.

SHEET TITLE: TABLE OF CONTENTS



NC.
TEL (916) 920-2020 Date: 06-10-2016
www.cyseng.com Page: 1 of 14

LTScale:6 Time:Jun10,2016-03:14pm Login:mayerhoferm Dimscale:1 BACTEC MGIT 320 OPM-0345-13\STRU\S1_TASK 03.dwq 8 OPMs\Task 03 Three L:\Jobs15\15117

BD Life Sciences - Diagnostic Systems BACTECTM MGITTM 320



GENERAL NOTES:

- 1. THIS OSHPD PRE-APPROVAL OF MANUFACTURER'S CERTIFICATION (OPM) IS BASED ON THE CBC 2013. THE DEMAND (DESIGN FORCES) FOR USE WITH THIS OPM SHALL BE BASED ON THE CBC 2013.
- 2. IT IS THE RESPONSIBILITY OF THE SEOR FOR A SITE SPECIFIC PROJECT TO VERIFY:
 - A. THE ADEQUACY OF THE NEW OR EXISTING STRUCTURE TO RESIST THE FORCES & WEIGHT SPECIFIED FOR EACH EQUIPMENT IN ADDITION TO ALL OTHER LOADS. PROVIDE & DESIGN SUPPLEMENTARY MEMBERS AS REQUIRED.
 - B. THAT THE FLOOR ANCHORS ARE LOCATED AT AN ADEQUATE DISTANCE FROM ANY SLAB EDGES OR OPENINGS.
 - C. THAT THE FLOOR ANCHORS ARE LOCATED AT AN ADEQUATE DISTANCE FROM ANY NEW OR EXISTING ANCHORS. THE SPACING SHOWN IN THE TEST LOADS TABLE ON <u>PAGE 3</u> IS THE REQUIRED MINIMUM SPACING OF THE GIVEN DIAMETER ANCHORS. THE REQUIRED SPACING FROM ANCHORS OF OTHER DIAMETERS & EMBEDMENTS MAY VARY & SHALL BE EVALUATED BY THE SEOR.
 - D. THAT THE INSTALLATION IS IN CONFORMANCE WITH THE CBC 2013 & WITH THE DETAILS SHOWN IN THIS PRE-APPROVAL.
 - E. THAT THE ACTUAL EQUIPMENT'S WEIGHT, CENTER OF GRAVITY (CG) LOCATION, ANCHOR LOCATIONS, ANCHOR DETAILS, & THE MATERIAL & GAGE OF THE EQUIPMENT WHERE ATTACHMENTS ARE MADE, AGREE WITH THE INFORMATION SHOWN ON THE PRE-APPROVAL DOCUMENTS.
- 3. EXPANSION ANCHORS INSTALLED IN NORMAL WEIGHT OR SAND-LIGHTWEIGHT CONCRETE SHALL BE CARBON STEEL HILTI KB-TZ EXPANSION ANCHORS COMPLYING WITH ESR-1917 REISSUED MAY 2015.

 A. INSTALLATION: INSTALL THE EXPANSION ANCHORS IN ACCORDANCE WITH THE REQUIREMENTS GIVEN IN
 - A. INSTALLATION: INSTALL THE EXPANSION ANCHORS IN ACCORDANCE WITH THE REQUIREMENTS GIVEN IN THE ICC EVALUATION REPORT FOR THE SPECIFIC ANCHOR & THE PARAMETERS GIVEN IN THE TABLE ON PAGE 3.
 - B. JOB TESTING: FOR VERIFYING SATISFACTORY INSTALLATION WORKMANSHIP, PERFORM JOB SITE TESTING IN ACCORDANCE WITH THE TEST LOAD TABLE PROVIDED IN THIS DOCUMENT. TEST 50% OF THE INSTALLED ANCHORS. THE TEST LOAD MAY BE APPLIED BY ANY METHOD THAT WILL EFFECTIVELY MEASURE THE TORQUE IN THE ANCHOR SUCH AS CALIBRATED TORQUE WRENCH METHOD. ALL TESTS SHALL BE CONDUCTED IN THE PRESENCE OF THE INSPECTOR OF RECORD (IOR). IF ANY ANCHOR FAILS THE TEST, TEST ALL ANCHORS. THE TEST SHALL BE PERFORMED 24 HOURS OR MORE AFTER INSTALLATION. TESTING MAY BE DONE PRIOR TO EQUIPMENT INSTALLATION. ALSO REFER TO CBC 1913A.7 "TESTS FOR POST—INSTALLED ANCHORS IN CONCRETE".
 - C. FAILURE/ACCEPTANCE CRITERIA: THE FOLLOWING CRITERIA APPLY FOR THE ACCEPTANCE OF INSTALLED ANCHORS:
 - TORQUE WRENCH METHOD: THE APPLICABLE TEST TORQUE MUST BE REACHED WITHIN THE FOLLOWING LIMITS:


WEDGE TYPE: ONE-HALF (1/2) TURN OF THE NUT.



SHEET TITLE: GENERAL NOTES

CYS STRUCTURAL ENGINEERS, INC.	No: 15117.03
2495 NATOMAS PARK DRIVE, SUITE 650 TEL (916) 920-2020 Date	: 06-10-2016
SACRAMENTO, CA 95833 www.cyseng.com Page	e: 2 of 14

GENERAL NOTES CONTINUED:



ANCHOR DIA (INCH)	INSTALLATION EMBED (INCH) hnom	EFFECTIVE EMBED (INCH) heff	HOLE DEPTH (INCH) ho	MIN CONC THICKNESS (INCH) h	MIN CONC EDGE DISTANCE (INCH)	MIN ANCHOR SPACING (INCH)	TEST TORQUE (FT-LBS)	CONDITION OF ANCHORAGE
3/8*	25/16	2	25/8	C C	OD.	2.50	25	CASE 1*
1/2	23/8	2	25/8	4	12 0,	5.00	40	CASE 2

- * STRUT & SUPPORT ANCHORS IN THE SOFFIT OF CONCRETE OVER METAL DECK
- BOLTS THROUGH CONCRETE ON METAL DECK: 345-13

 A. BOLTS SHALL BE TORQUED BY ¾ TURN OF THE NUT AFTER SNUG TIGHT CONDITION IS ACHIEVED, UND. THE SNUG TIGHT CONDITION IS DEFINED AS THE TIGHTNESS REQUIRED TO BRING THE CONNECTED PLIES INTO FIRM CONTACT ffrey Y. Kikum

 - THRU-BOLT HOLES SHALL BE \mathcal{H}_6 " LARGER THAN BOLT SIZE (HOLE SIZE = BOLT SIZE + \mathcal{H}_6 "). THRU-BOLTS IN CONCRETE SHALL RECEIVE SPECIAL INSPECTION & TESTING IN ACCORDANCE WITH REQUIREMENTS FOR POST-INSTALLED FANCHORS! THRU-BOLTS WITH STEEL TO STEEL CONNECTION IN TENSION DO NOT REQUIRE TESTING.

SHEET TITLE: GENERAL NOTES (CONTINUED)

CYS STRUCTURAL ENGINEERS, INC.

2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833


TEL (916) 920-2020

Job No: 15117.03 06-10-2016 Date: 3 of 14 www.cyseng.com Page.



GENERAL NOTES CONTINUED:

TWO (2) CASES OF ATTACHMENT ARE SPECIFIED & PRESENTED IN THIS PRE-APPROVAL:



CASE 1: ATTACHMENT DETAILS LOCATED AT UPPER FLOORS ABOVE THE BASE OF A BLDG ($z/h \le 1$). THE FLOORS ARE ASSUMED TO BE BUILT OF A MIN 31/4" SLWC TOPPING OVER MIN 20 GA MTL DECK (f'c = 3000 PSI, MIN).

CASE 2: ATTACHMENT DETAILS LOCATED AT OR BELOW THE BASE OF A BLDG (z/h=0). THE FLOORS ARE ASSUMED TO BE BUILT OF A MIN 4" NWC SLAB (f'c = 3000 PSI, MIN).

THIS PRE-APPROVAL MAY BE USED AT ANY GEOGRAPHICAL LOCATION IN THE STATE OF CALIFORNIA. WHERE S_{DS} IS LESS THAN OR EQUAL TO 2.50.

SYSTEM OVERVIEW & DESIGN CRITERIA:

- A SINGLE INSTALLATION CONFIGURATION OF THE MGIT 320 INSTRUMENT IS COVERED BY THIS OPM, THE STAND CONFIGURATION.
- 2. SUPPORT & ATTACHMENT DESIGN IS PER 2013 CBC AT LRFD LEVEL FORCES.

OTHER MECHANICAL OR ELECTRICAL COMPONENTS PER TABLE 13.6-1 OF ASCE 7-10 SUPPLEMENT #1:

$$a_p = 1.0$$
 $R_p = 1.5$ DATI = 1.514/2010 = 1.5 (FOR CONC ANCHORS ONLY)


W_P AS NOTED ON DRAWING ON PG 6

UPPER FLOORS ABOVE THE BASE OF BUILDING, z/h = 1 $S_{DS} = 2.50$ $F_{D} = 3.00$ W_{D} CASE 1:

FLOORS AT OR BELOW THE BASE OF BUILDING, z/h = 0 CASE 2: $S_{0S} = 2.50$ $F_{\rm p} = 1.13 \, W_{\rm p}$

LOAD COMBINATIONS

 $(0.9 - 0.2 S_{DS}) D - \Omega_0 F_p$ (FOR MAXIMUM TENSION) $(1.2 + 0.2 S_{DS}) D + \Omega_0 F_D$ (FOR MAXIMUM COMPRESSION)



SHEET TITLE: GENERAL NOTES (CONTINUED)

SYSTEM OVERVIEW & DESIGN CRITERIA



CYS STRUCTURAL ENGINEERS, INC.

2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833

(916) 920-2020 Date:

Job No: 15117.03 06-10-2016 4 of 14 www.cyseng.com Page:



ABBREVIATIONS:


0	AT ANOLIOD BOLT	f'c	MINIMUM ULTIMATE COMPRESSIVE	OPG	OPENING
AB ABV	ANCHOR BOLT ABOVE	EL D	STRENGTH OF CONCRETE	PERP	PERPENDICULAR PACE(S)
ADJ	ADJACENT	FLR FT (')	FLOOR FOOT/FEET	PG(S) PL	PAGE(S) PLATE
ALUM	ALUMINUM	F _D	HORIZONTAL SEISMIC FORCE	PSI	POUNDS PER SQUARE INCH
ASTM	AMERICAN SOCIETY FOR	' P	PER ASCE 7–10 SEISMIC	RECT	RECTANGULAR
	TESTING & MATERIALS		FORCE REQUIREMENTS	SEOR	STRUCTURAL ENGINEER OF
BLDG	BUILDING	F_{y}	SPECIFIED MINIMUM YIELD	020	RECORD
BLW	BELOW	•	STRESS OF STEEL	SLWC	SAND-LIGHTWEIGHT CONCRETE
BOTT	BOTTOM PHILIPING CODE	GA	GAUGE	SS	STAINLESS STEEL
CBC CG	CALIFORNIA BUILDING CODE CENTER OF GRAVITY	ICC	INTERNATIONAL CODE	STL	STEEL
Ę	CENTERLINE	INI /"\	COUNCIL	THRD	THREAD OR THREADED
CONC	CONCRETE	IN (") KSI	INCH KIPS PER SQUARE INCH	Tu	ANCHORAGE TENSION REACTION DUE TO SEISMIC FORCE
COORD	COORDINATE	l VOI	ANGLE	TYP	TYPICAL
CRS	COLD ROLLED STEEL	LBS	POUNDS	T&B	TOP & BOTTOM
DIA (ø)	DIAMETER	LRFD	LOAD AND RESISTANCE	Vu	ANCHORAGE SHEAR REACTION
(E)	EXISTING CONDITION		FACTOR DESIGN E		DUE TO SEISMIC FORCE
EA	EACH	MAX	MAXIMUM	W/	WITH
EE	EACH END	MFR	MANUFACTURER	Wp	OPERATING WEIGHT
ELEV	ELEVATION	MIN	MINIMUM DO DO METAD	WT	WEIGHT
EQ EQUIP	EQUIPMENT	NO. (#)			
ES	EACH SIDE		NORMAL4-WEIGHT CONCRETE	Z	
		11110	THOMMAL MEIOTIP CONTONLIL		

BY: Jeffrey Y. Kikumoto


DATE: 07/14/2016

PROFESS IONAL TO SHARE TO SHARE OF CALIFORNIA TO STAUCTURAL TO STAUCTURA

SHEET TITLE: ABBREVIATIONS







SHEET TITLE: STAND CONFIGURATION **ELEVATIONS**

CYS STRUCTURAL ENGINEERS, INC.

2495 NATOMAS PARK DRIVE, SUITE 650

(916) 920-2020 Date: TEL www.cyseng.com Page.


Job No: 15117.03 06-10-2016 6 of 14

SACRAMENTO, CA 95833



	М	AX LRF	D FORG	CES AT	LEVELIN	NG LEG	S ¹ (LBS	S)
FRONT BRACKET REAR BRAC							RACKET	
	LEC	G 1	LEG	9 2	LEG	3	LEC	4
	Tu	Vux	Tu	Vux	Tu	Vux	Tu	Vux
CASE 1	180	260	180	260	180	260	180	260
CASE 2	50	100	50	100	50	100	50	100

. ECCENTRICITY & PRYING ACTION MUST BE CONSIDERED BASED ON THE BRACKET CONFIGURATION.



SHEET TITLE: ATTACHMENT PLAN UNIT TO STAND

CYS STRUCTURAL ENGINEERS, INC.

2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833 TEL (916) 920-2020 Date: www.cyseng.com Page:

Job No: 15117.03 Date: 06-10-2016 Page: 7 of 14

34.55"

32.05"

ASTM A36 STL

P WASHER 3/8×2/2×0'-4"

(

0

W/ STD SIZE HOLE

AT EA THRU-BOLT

0

FRONT


PLAN VIEW



	MAX LRFD FORCES AT LEVELING LEGS ¹ (LBS)							
		FRONT E	BRACKET			REAR B	RACKET	
	LEG 1 LEG 2			LEC	3	LEG 4		
	Tu	$\Omega_{o}V_{ux}$	Tu	ΩοVux	Tu	ΩοVux	Tu	ΩοVux
CASE 1	730	160	670	230	0	490	1570	750
	Ω _o Tu	$\Omega_{o}V_{ux}$	Ω _o Tu	ΩοVux	Ω _o Tu	ΩοVux	Ω _o Tu	ΩοVux
CASE 2	380	60	350	90	860	200	0	270

1.25"

- **ECCENTRICITY & PRYING ACTION** MUST BE CONSIDERED BASED ON THE BRACKET CONFIGURATION.
- OVERSTRENGTH FACTOR $(\Omega_{\rm o})$ INCLUDED WHERE NOTED.



NOTE:

DO NOT USE INTERIOR HOLES FOR THRU-BOLTED CONDITION AT SUSPENDED FLR, TYP EA BRACKET

FRONT MOUNT HEAVY ANCHOR PER PG 9

STAND OUTLINE

½"ø LEVELING LEG (TYP OF 4)

¾"ø HEX CAP SCREW ₩∠ WASHER & LOCKWASHER SUPPLIED BY MFR (TYP OF 4)

NOTE: SEE PG 6 FOR CG LOCATION & WT.

ATE OF CALL

SHEET TITLE: ATTACHMENT PLAN STAND TO FLOOR

0

0

0

CYS STRUCTURAL ENGINEERS, INC.

2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833

(916) 920-2020 Date: TEL www.cyseng.com Page.

Job No: 15117.03 06-10-2016 8 of 14

BACTEC MGIT 320 OPM-0345-13\STRU\S1_TASK 03.dwg Time:Jun10,2016-03:15pm Login:mayerhoferm Dimecale:1 LTScale:6

8


Three OPMs\Task 03


:\Jobs15\15117

33.38" 88


30.

1.25"









SHEET TITLE: LEVELING LEG DETAILS


	CYS STRUCTURAL ENGINEERS,	INC.		Job No:	15117.03
	2495 NATOMAS PARK DRIVE, SUITE 650	TEL	(916) 920-2020	Date:	06-10-2016
l	SACRAMENTO, CA 95833		www.cyseng.com	Page:	11 of 14

Three OPMs\Task 03 BD BACTEC MGIT 320 OPM-0345-13\STRU\S1_TASK 03.dwg Time:Jun10,2016-03:17pm Login:mayerhoferm Dimscale:1 LTScale:6

L:\Jobs15\15117

BD Life Sciences - Diagnostic Systems BACTECTM MGITTM 320





CYS STRUCTURAL ENGINEERS, INC.


2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833

(916) 920-2020 Date: TEL www.cyseng.com Page.

Job No: 15117.03 06-10-2016 12 of 14

BD Life Sciences - Diagnostic Systems BACTECTM MGITTM 320






CYS STRUCTURAL ENGINEERS, INC. Job No: 15117.03 (916) 920-2020 Date: 06-10-2016 TEL 2495 NATOMAS PARK DRIVE, SUITE 650 13 of 14 www.cyseng.com | Page: SACRAMENTO, CA 95833



BACTEC MGIT 320 OPM-0345-13\STRU\S1_TASK 03.dwq Time:Jun10,2016-03:18pm Login:mayerhoferm Dimscale:1 LTScale:6 8 Three OPMs\Task 03



BUILDING



STAND TO 4" SLAB ON GRADE (CASE 2)

ı

L:\Jobs15\15117

CYS STRUCTURAL ENGINEERS, ÍNC. 2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833

(916) 920-2020 Date: TEL

Job No: 15117.03 06-10-2016 www.cyseng.com Page: 14 of 14