



OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT
FACILITIES DEVELOPMENT DIVISION

APPLICATION FOR OSHPD PREAPPROVAL OF
MANUFACTURER'S CERTIFICATION (OPM)

OFFICE USE ONLY

APPLICATION #: OPM-0492

OSHPD Preapproval of Manufacturer's Certification (OPM)

Type: ☒ New ☐ Renewal/Update

Manufacturer Information

Manufacturer: Becton, Dickinson and Company

Manufacturer's Technical Representative: Dustin Diemert

Mailing Address: 7 Loveton Circle, Sparks, MD 21152

Telephone: (410) 316-4862

Email: dustin_diemert@bd.com

Product Information

Product Name: BD CORTM SYSTEM

Product Type: High-Throughput Molecular Platform

Product Model Number: COR GX, COR MX, COR PX

General Description: An automated diagnostic instrument for use in clinical laboratories that supports a menu of clinically differentiated assays for women's health, sexually transmitted infections, and gastrointestinal (GI) applications.

Applicant Information

Applicant Company Name: Becton, Dickinson and Company

Contact Person: Dustin Diemert

Mailing Address: 7 Loveton Circle, Sparks, MD 21152

Telephone: (410) 316-4862

Email: dustin_diemert@bd.com

Title: Senior Engineer - Mechanical

"Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"

STATE OF CALIFORNIA – HEALTH AND HUMAN SERVICES AGENCY

OSHPD



**OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT
FACILITIES DEVELOPMENT DIVISION**

Registered Design Professional Preparing Engineering Recommendations

Company Name: CYS STRUCTURAL ENGINEERS, INC.

Name: Dieter Siebald

California License Number: S4346

Mailing Address: 2495 Natomas Park Drive, Suite 650, , Sacramento, CA 95833

Telephone: (916) 920-2020

Email: dieters@cyseng.com

OSHDP Special Seismic Certification Preapproval (OSP)

☐ Special Seismic Certification is preapproved under OSP

OSP Number: _____

Certification Method

Testing in accordance with: ☐ ICC-ES AC156 ☐ FM 1950-16

☐ Other(s) (Please Specify): _____

*Use of criteria other than those adopted by the California Building Standards Code, 2019 (CBSC 2019) for component supports and attachments are not permitted. For distribution system, interior partition wall, and suspended ceiling seismic bracings, test criteria other than those adopted in the CBSC 2019 may be used when approved by OSHDP prior to testing.

☒ Analysis

☐ Experience Data

☐ Combination of Testing, Analysis, and/or Experience Data (Please Specify): _____

OSHDP Approval

Date: 7/20/2020

Name: William Staehlin

Title: Senior Structural Engineer

Condition of Approval (if applicable): _____




C:\Users\camachom\appdata\local\temp\AcPublish_28780\51.dwg Time:Jun25,2020--08:42am Login:camachom Dimstyle:1 LTScale:6

TABLE OF CONTENTS
OPM-0492-19

	PAGE
GENERAL NOTES	2
ABBREVIATIONS, SYSTEM OVERVIEW & DESIGN CRITERIA	3
COR CONFIGURATION PLAN & ELEVATIONS	4
SEISMIC BRACKET ASSEMBLIES	5
SEISMIC BRACKET DETAILS	6
ATTACHMENT DETAILS	
CONCRETE FILL OVER METAL DECK (CASE 1)	7
6" CONCRETE SLAB ON GRADE (CASE 2)	8
4" CONCRETE SLAB ON GRADE (CASE 3)	8

- NOTES:**
- 1. THESE DRAWINGS ARE PREPARED FOR BD LIFE SCIENCES – DIAGNOSTIC SYSTEMS, SPARKS, MARYLAND.
 - 2. THE CONTRACTOR & INSPECTOR OF RECORD SHALL OBTAIN A COPY OF THIS PRE-APPROVAL FROM THE OFFICE OF STATEWIDE HEALTH PLANNING & DEVELOPMENT (OSHDP) PRE-APPROVAL PROGRAMS WEBSITE.
 - 3. THIS PRE-APPROVAL COVERS THE SUPPORTS & ATTACHMENTS OF THE EQUIPMENT TO THE SUPPORTING STRUCTURE. THE EQUIPMENT & ATTACHMENT HARDWARE ARE SUPPLIED BY THE MANUFACTURER. THE EXPANSION ANCHORS, THRU-BOLTS & STRUT PLATES SHOWN IN THIS OPM SHALL BE SUPPLIED & INSTALLED BY THE CONTRACTOR.




SHEET TITLE: TABLE OF CONTENTS				Rev	Description	Date	Job No: 20041
<div><div>BD Life Sciences - Diagnostic Systems BD COR™</div></div> <div><div>CYS STRUCTURAL ENGINEERS, INC. 2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833</div></div> <div><div>TEL (916) 920-2020 www.cyseng.com</div></div>							Date: 6-25-2020
							By: MTC
							Page: 1 of 8

GENERAL NOTES:


1. THIS OSHPD PRE-APPROVAL OF MANUFACTURER'S CERTIFICATION (OPM) IS BASED ON THE CBC 2019. THE DEMAND (DESIGN FORCES) FOR USE W/ THIS OPM SHALL BE BASED ON THE CBC 2019.
2. IT IS THE RESPONSIBILITY OF THE SEOR FOR A SITE SPECIFIC PROJECT TO VERIFY:
- A. THE ADEQUACY OF THE NEW OR EXISTING STRUCTURE TO RESIST THE FORCES & WT SPECIFIED FOR EA EQUIP IN ADDITION TO ALL OTHER LOADS. PROVIDE & DESIGN SUPPLEMENTARY MEMBERS AS REQ.
 - B. THAT THE FLR ANCHORS ARE LOCATED AT AN ADEQUATE DISTANCE FROM ANY SLAB EDGES OR OPGS.
 - C. THAT THE FLR ANCHORS ARE LOCATED AT AN ADEQUATE DISTANCE FROM ANY NEW OR EXISTING ANCHORS. THE SPCG SHOWN IN THE TEST LOADS TABLE ON PG 2 IS THE REQ MIN SPCG OF THE GIVEN DIA ANCHORS. THE REQ SPCG FROM ANCHORS OF OTHER DIAMETERS & EMBEDMENTS MAY VARY & SHALL BE EVALUATED BY THE SEOR.
 - D. THAT THE INSTALLATION IS IN CONFORMANCE W/ THE CBC 2019 & W/ THE DETAILS SHOWN IN THIS PRE-APPROVAL.
 - E. THAT THE ACTUAL EQUIP'S WT, CENTER OF GRAVITY (CG) LOCATION, ANCHOR LOCATIONS, ANCHOR DETAILS, & THE MATERIAL & GAGE OF THE EQUIP WHERE ATTACHMENTS ARE MADE, AGREE W/ THE INFO SHOWN ON THE PRE-APPROVAL DOCUMENTS.
3. EXPANSION ANCHORS INSTALLED IN NWC OR SLWC SHALL BE CARBON STEEL HILTI KB-TZ EXPANSION ANCHORS COMPLYING W/ ESR-1917 REISSUED JANURAY 2020. PROVIDE FULL THREAD ENGAGEMENT OF NUT & WASHER.
- A. INSTALLATION: INSTALL THE EXPANSION ANCHORS IN ACCORDANCE W/ THE REQUIREMENTS GIVEN IN THE ICC EVALUATION REPORT FOR THE SPECIFIC ANCHOR & THE PARAMETERS GIVEN IN THE TABLE ON PG 2.
 - B. JOB TESTING: FOR VERIFYING SATISFACTORY INSTALLATION WORKMANSHIP, PERFORM JOBSITE TESTING IN ACCORDANCE W/ THE TEST LOAD TABLE PROVIDED IN THIS DOCUMENT. TEST 50% OF THE INSTALLED ANCHORS. ALL TESTS SHALL BE CONDUCTED IN THE PRESENCE OF THE SPECIAL INSPECTOR. TESTING & SPECIAL INSPECTION OF EXPANSION ANCHORS SHALL BE PERFORMED BY THE FACILITY OWNER PER CBC 1704A & 1910A.5 & CAC 7-149. ALL REPORTS SHALL BE SENT TO THE INSPECTOR OF RECORD, OWNER & ARCHITECT OR ENGINEER IN RESPONSIBLE CHARGE. IF ANY ANCHOR FAILS THE TEST, TEST ALL ANCHORS. THE TEST SHALL BE PERFORMED 24 HOURS OR MORE AFTER INSTALLATION. TESTING MAY BE DONE PRIOR TO EQUIP INSTALLATION, HOWEVER, THE NUT SHALL BE RETORQUED TO INSTALLATION TORQUE AFTER EQUIP INSTALL. ALSO, REFER TO 2019 CBC 1910A.5 "TESTS FOR POST-INSTALLED ANCHORS IN CONCRETE".
 - C. FAILURE/ACCEPTANCE CRITERIA: THE FOLLOWING CRITERIA APPLY FOR THE ACCEPTANCE OF INSTALLED ANCHORS:

- **TORQUE WRENCH METHOD:** THE APPLICABLE TEST TORQUE MUST BE REACHED WITHIN THE FOLLOWING LIMITS:
WEDGE TYPE: ONE-HALF (1/2) TURN OF THE NUT.



CONDITION OF ANCHORAGE	ANCHOR DIA (INCH)	INSTALLATION EMBED (INCH) hnom	EFFECTIVE EMBED (INCH) hef	HOLE DEPTH (INCH) ho	MIN CONC THK (INCH) h	MIN CONC EDGE DISTANCE (INCH)	MIN ANCHOR SPCG (INCH)	TEST TORQUE (FT-LBS)	MIN SPCG TO (E) ANCHORS (INCH)
CASE 1 STRUT PL	3/8	2 5/16	2	2 5/8	3/4	12	6 3/4	25	5
CASE 2	1/2	3 5/8	3 1/4	3 3/8	6	12	4	40	5
CASE 3	1/2	2 3/8	2	2 5/8	4	12	4	40	5

- D. AVOID DAMAGING (E) STL REINF IN CONC SLAB WHEN INSTALLING CONC EXPANSION ANCHORS.
 - E. PROVIDE FOR FULL THRD ENGAGEMENT OF NUTS & WASHER.
4. BOLTS THRU CONC ON MTL DECK:
- A. BOLTS SHALL BE TORQUED BY 3/4 TURN OF THE NUT AFTER SNUG TIGHT CONDITION IS ACHIEVED, UNO. THE SNUG TIGHT CONDITION IS DEFINED AS THE TIGHTNESS REQ TO BRING THE CONNECTED PLIES INTO FIRM CONTACT.
 - B. THRU-BOLT HOLES SHALL BE 1/16" LARGER THAN BOLT SIZE (HOLE SIZE = BOLT SIZE + 1/16").
 - C. THRU-BOLTS IN CONC SHALL RECEIVE SPECIAL INSPECTION & TESTING IN ACCORDANCE W/ REQUIREMENTS FOR POST-INSTALLED ANCHORS. THRU-BOLTS W/ STL TO STL CONNECTION IN TENSION DO NOT REQUIRE TESTING.
5. THREE (3) CASES OF ATTACHMENT ARE SPECIFIED & PRESENTED IN THIS PRE-APPROVAL:



CASE 1: ATTACHMENT DETAILS LOCATED AT UPPER FLRS ABV THE BASE OF A BLDG. THE FLRS ARE ASSUMED TO BE BUILT OF A MIN 3 1/4" SLWC TOPPING OVER 3" DEEP MIN 20 GA MTL DECK (f'c = 3000 PSI, MIN). $z/h \leq 1$ & $S_{ps} \leq 2.50$

CASE 2: ATTACHMENT DETAILS LOCATED AT OR BLW THE BASE OF A BLDG. THE FLRS ARE ASSUMED TO BE BUILT OF A MIN 6" NWC SLAB (f'c = 3000 PSI, MIN). $z/h = 0$ & $S_{ps} \leq 2.50$

CASE 3: ATTACHMENT DETAILS LOCATED AT OR BLW THE BASE OF A BLDG. THE FLRS ARE ASSUMED TO BE BUILT OF A MIN 4" NWC SLAB (f'c = 3000 PSI, MIN). $z/h = 0$ & $S_{ps} \leq 1.8$

6. THIS PRE-APPROVAL MAY BE USED ONLY AT GEOGRAPHICAL LOCATIONS IN THE STATE OF CALIFORNIA WHERE S_{ps} & z/h IS LESS THAN OR EQ TO THE VALUES NOTED ABV.



SHEET TITLE: GENERAL NOTES				Rev	Description	Date	Job No: 20041
<div>BD Life Sciences - Diagnostic Systems BD COR™</div>				<div>CYS STRUCTURAL ENGINEERS, INC. 2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833</div>		TEL (916) 920-2020 www.cyseng.com	Date: 6-25-2020
							By: MTC
							Page: 2 of 8

C:\Users\camachom\appdata\local\temp\AcPublish_28780\51.dwg Time:Jun25,2020-08:42am Login:camachom Dimstyle:1 LTScale:6

ABBREVIATIONS:

@	AT	f'c	MINIMUM ULTIMATE COMPRESSIVE STRENGTH OF CONCRETE
∠	ANGLE		
AB	ANCHOR BOLT	FLG	FLANGE
ABV	ABOVE	FLR	FLOOR
ASCE	AMERICAN SOCIETY OF CIVIL ENGINEERS	FT (')	FOOT/FEET
		F _p	HORIZONTAL SEISMIC FORCE PER ASCE 7-10 SEISMIC FORCE REQUIREMENTS
ADJ	ADJACENT	F _y	SPECIFIED MINIMUM YIELD STRESS OF STEEL
AISI	AMERICAN IRON & STEEL INSTITUTE		
ALUM	ALUMINUM	GA	GAUGE
ASTM	AMERICAN SOCIETY FOR TESTING & MATERIALS	GR	GRADE
BLDG	BUILDING	HT	HEIGHT
BLW	BELOW	ICC	INTERNATIONAL CODE COUNCIL
BOTT	BOTTOM		
CBC	CALIFORNIA BUILDING CODE	IN (")	INCH
CG	CENTER OF GRAVITY	INFO	INFORMATION
CL	CENTERLINE	KSI	KIPS PER SQUARE INCH
CONC	CONCRETE	LBS	POUNDS
CONT	CONTINUOUS	LRFD	LOAD AND RESISTANCE FACTOR DESIGN
COORD	COORDINATE		
CRS	COLD-ROLLED STEEL	MAX	MAXIMUM
DBL	DOUBLE	MFR	MANUFACTURER
DEG	DEGREE	MIN	MINIMUM
DIA (ø)	DIAMETER	mm	MILLIMETER
DTL	DETAIL	MTL	METAL
(E)	EXISTING CONDITION	NO. (#)	NUMBER OR POUNDS
EA	EACH	NWC	NORMAL WEIGHT CONCRETE
EE	EACH END		
ELEV	ELEVATION		
EQ	EQUAL		
EQUIP	EQUIPMENT		
ES	EACH SIDE		
EXTR	EXTERIOR		

OP	OPERATING
OPG	OPENING
OPM	OSPHD PRE-APPROVAL OF MANUFACTURER'S CERTIFICATION
	OFFICE OF STATEWIDE HEALTH PLANNING & DEVELOPMENT
OSPHD	PERPENDICULAR
PERP	PAGE
PG	PLATE
ℓ	POUNDS PER SQUARE INCH
PSI	REQUIRED
REQ	STRUCTURAL ENGINEER OF RECORD
SEOR	SAND-LIGHTWEIGHT CONCRETE
	SPACING
SLWC	SPECIFICATION
SPCG	STAINLESS STEEL
SPEC	STEEL
SS	THICK/THICKNESS
STL	THREAD OR THREADED
THK	TOP OF
THRD	ANCHORAGE TENSION REACTION DUE TO SEISMIC FORCE
T.O.	TYPICAL
Tu	TOP & BOTTOM
	UNLESS NOTED OTHERWISE
TYP	VERTICAL
T&B	ANCHORAGE SHEAR REACTION DUE TO SEISMIC FORCE
UNO	WITH
VERT	OPERATING WEIGHT
Vu	WEIGHT
W/	
W _p	
WT	

SYSTEM OVERVIEW & DESIGN CRITERIA:

1. THE COR UNIT INCLUDES THREE (3) MODULES THAT CAN BE INTERCONNECTED (COR PX, COR GX & COR MX). INSTALLATION CONFIGURATIONS ARE SHOWN IN TABLE BLW:

CONFIGURATION	LEFT MODULE	CENTER MODULE	RIGHT MODULE
1	COR GX	COR PX	COR GX
2	COR GX	COR PX	—
3	—	COR PX	COR GX
4	COR MX	COR PX	—
5	—	COR PX	COR MX
6	COR MX	COR PX	COR MX
7	COR GX	COR PX	COR MX
8	COR MX	COR PX	COR GX

2. SUPPORT & ATTACHMENT DESIGN IS PER 2019 CBC AT LRFD LEVEL FORCES.

OTHER MECHANICAL OR ELECTRICAL COMPONENTS PER TABLE 13.6-1 OF ASCE 7-16 SUPPLEMENT #1:

$\alpha_p = 1.0$ $R_p = 1.5$ $I_p = 1.5$ $\Omega_0 = 1.5$ (FOR CONC ANCHORS ONLY)

$W_p = 1573\#$ MAX PER MODULE

UPPER FLRS ABV THE BASE OF BLDG, $z/h \leq 1$

CASE 1: $S_{Ds} = 2.50$ $F_p = 3.00 W_p$

FLRS AT OR BLW THE BASE OF BLDG, $z/h \leq 0$


CASE 2: $S_{Ds} = 2.50$ $F_p = 1.13 W_p$

CASE 3: $S_{Ds} = 1.8$ $F_p = 0.81 W_p$

LOAD COMBINATIONS

$(0.9 - 0.2 S_{Ds}) D - \Omega_0 F_p$ (FOR MAX TENSION)

$(1.2 + 0.2 S_{Ds}) D + \Omega_0 F_p$ (FOR MAX COMPRESSION)



SHEET TITLE: ABBREVIATIONS, SYSTEM OVERVIEW & DESIGN CRITERIA



BD Life Sciences - Diagnostic Systems
BD COR™




CYS STRUCTURAL ENGINEERS, INC.

2495 NATOMAS PARK DRIVE, SUITE 650
SACRAMENTO, CA 95833

TEL (916) 920-2020
www.cyseng.com


Rev	Description	Date	Job No: 20041
			Date: 6-25-2020
			By: MTC
			Page: 3 of 8

C:\Users\camachom\appdata\local\temp\AcPublish_28780\51.dwg Time:Jun25,2020-08:42am Login:camachom Dimscale:1 LScale:6



CG LOCATION & MAX OP WT				
INSTRUMENT	WT (LB)	X (mm)	Y (mm)	Z (mm)
CENTER MODULE	1405	30 1/2" (776)	12 5/8" (320)	23 7/8" (607)
LEFT MODULE COR GX	1573	30 3/8" (772)	13 5/8" (347)	28 1/8" (715)
RIGHT MODULE COR GX	1573	30 3/8" (772)	13 5/8" (347)	28 1/8" (715)
LEFT MODULE COR MX	1416	31 7/8" (810)	13 5/8" (346)	26 15/16" (684)
RIGHT MODULE COR MX	1416	31 7/8" (810)	13 5/8" (346)	26 15/16" (684)

NOTE:
ONE "COR GX" OR "COR MX" MODULE SHALL OCCUR ON EITHER THE RIGHT, LEFT OR BOTH SIDES OF THE "COR PX" CENTER MODULE.



SHEET TITLE: COR CONFIGURATION PLAN & ELEVATIONS		Rev	Description	Date	Job No: 20041
 BD Life Sciences - Diagnostic Systems BD COR™					Date: 6-25-2020
					By: MTC
					Page: 4 of 8




BD Life Sciences - Diagnostic Systems
BD COR™




CYS STRUCTURAL ENGINEERS, INC.
2495 NATOMAS PARK DRIVE, SUITE 650
SACRAMENTO, CA 95833

TEL (916) 920-2020
www.cyseng.com




NOTE:
LEVELING PAD BASE IS NOT
CONNECTED TO THE FLR BRACKET.


SECTION D-D




SECTION A-A



SECTION B-B




SECTION C-C




SHEET TITLE: SEISMIC BRACKET ASSEMBLIES

	BD Life Sciences - Diagnostic Systems BD COR™		 CYS STRUCTURAL ENGINEERS, INC. 2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833 TEL (916) 920-2020 www.cyseng.com		Rev	Description	Date	Job No: 20041
								Date: 6-25-2020
								By: MTC
								Page: 5 of 8

L:\Jobs\20\20041 BD - Amendment to OPM-0492-13\ACAD\STRU\S1.dwg Time:Jun25,2020-08:57am Login:camachom DimScale:1 LTScale:6




PLAN VIEW




ELEV VIEW

1
—
1 1/2" = 1'-0"


FLOOR BRACKET
DETAIL



PLAN VIEW




ELEV VIEW




SIDE ELEV

2
—
1 1/2" = 1'-0"

INTERCONNECTION BRACKET
DETAIL




PLAN VIEW



SECTION A-A

3
—
6" = 1'-0"


LEVELING PAD BASE
DETAIL




SHEET TITLE: SEISMIC BRACKET DETAILS

	BD Life Sciences - Diagnostic Systems BD COR™		CYS STRUCTURAL ENGINEERS, INC. 2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833 TEL (916) 920-2020 www.cyseng.com	Rev	Description	Date	Job No: 20041
							Date: 6-25-2020
							By: MTC
							Page: 6 of 8


c:\Users\camachom\appdata\local\temp\Jun25,2020-08:42am Login:camachom Dimscale:1 L1Scale:6



SUSPENDED FLR W/ THRU-BOLTS




ANCHOR OPTIONS



PLAN VIEW

BRACKET(S) PERP TO FLUTES
AT EXTR BRACKETS




STRUT \bar{r} SUPPORT ANCHORS:
HILTI KB-TZ $\frac{3}{8}$ " ϕ x2" EMBED
(ICC ESR-1917) 1 EE OF \bar{r} MIN
(TORQUE TEST 25 FT-LBS)

NOTE:
USE \bar{r} IF THERE IS ONLY ONE BRACKET ABV.

PLAN VIEW

BRACKET(S) PARALLEL TO FLUTES

EXTEND STRUT TO NEXT ADJ LOW FLUTE
IF THRU-BOLTS ARE LESS THAN 2" FROM
STRUT \bar{r} SUPPORT ANCHORS



STRUT \bar{r} SUPPORT ANCHORS:
HILTI KB-TZ $\frac{3}{8}$ " ϕ x2" EMBED
(ICC ESR-1917) 2 EE OF \bar{r} MIN
(TORQUE TEST 25 FT-LBS)

THRU-BOLT PER
ELEV, TYP

BOTT OF MTL DECK
(LOW FLUTES SHADED
FOR CLARITY)

PLAN VIEW

BRACKET(S) PERP TO FLUTES
AT INTR BRACKETS

$S_{ps} \leq 2.5$ CASE 1 $z/h < 1$ W/O Ω_o	MAX LRFD FORCES AT EA ANCHOR	
	T_u	V_u
	2359#	860#
CASE 1 $z/h < 1$ W/ Ω_o	N/A	1290#

OVERSTRENGTH FACTOR (Ω_o) INCLUDED
ONLY IN SHEAR VALUE.



SHEET TITLE: ATTACHMENT DETAILS
CONCRETE FILL OVER METAL DECK (CASE 1)





BD Life Sciences - Diagnostic Systems
BD CORTM



CYS STRUCTURAL ENGINEERS, INC.
2495 NATOMAS PARK DRIVE, SUITE 650
SACRAMENTO, CA 95833
TEL (916) 920-2020
www.cyseng.com

Rev	Description	Date	Job No:	20041
			Date:	6-25-2020
			By:	MTC
			Page:	7 of 8

z/h = 0	MAX LRFD FORCES AT EA ANCHOR	
	T _u	V _u
CASE 2 S _{ps} ≤ 2.5	1301#	442#
CASE 3 S _{ps} ≤ 1.8	895#	318#



SHEET TITLE: ATTACHMENT DETAILS				Rev	Description	Date	Job No: 20041
6" CONCRETE SLAB ON GRADE (CASE 2) & 4" CONCRETE SLAB ON GRADE (CASE 3)							Date: 6-25-2020
 BD Life Sciences - Diagnostic Systems BD COR™	 CYS STRUCTURAL ENGINEERS, INC. 2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833		TEL (916) 920-2020 www.cyseng.com				By: MTC
							Page: 8 of 8