

OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT **FACILITIES DEVELOPMENT DIVISION**

APPLICATION FOR OSHPD PREAPPROVAL OF	OFFICE USE ONLY				
MANUFACTURER'S CERTIFICATION (OPM)	APPLICATION #: OPM-0208				
OSHPD Preapproval of Manufacturer's Certification (OPM)					
Type: New X Renewal/Update					
Manufacturer Information					
Manufacturer: BD					
Manufacturer's Technical Representative: Jared Zamaloff					
Mailing Address: 10020 Pacific Mesa Blvd., San Diego, CA 92121					
Telephone: (818) 876-4287 Email: jared.zamaloff@caref	usion.com				
FOR CODE CON					
Product Information OSHPD					
Product Name: PYXIS DUOSTATION	Z				
Product Type: Other Mechanical & Electrical Equipment	CH				
Product Model Number: Pyxis DuoStation By: Jeffrey Kikumoto					
General Description: Pyxis DuoStation system allows clinicians to safely and accumulation system system allows clinicians to safely and accumulation system system allows clinicians to safely and accumulation system syst	curately dispense both medications and medical				
	2				
Applicant Information	>				
Applicant Company Name: EASE LLC.					

Contact Person: Tiffany Tonn

Mailing Address: 1515 FAIRVIEW AVE, STE 205 MISSOULA, MT 59801

Telephone: (406) 541-3273 Email: tiffany@easeco.com

Title:

"Access to Safe. Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"

STATE OF CALIFORNIA- HEALTH AND HUMAN SERVICES AGENCY







OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT FACILITIES DEVELOPMENT DIVISION

Registered Design Professonal Preparing Engineering Recommendations							
Company Name: EASE LLC							
Name: Jonathan Roberson California License Number: S4197							
Mailing Address: 5877 Pine Ave., Suite 210, Chino Hills, CA 91709							
Felephone: (909) 606-7622 Email: jon@easeco.com							
OSHPD Special Seismic Certification Preapproval (OSP)							
Special Seismic Certification is preapproved under OSP OSP Number:							
OR GODE							
Certification Method							
Testing in accordance with:							
Other(s) (Please Specify):							
*Use of criteria other than those adopted by the California Building Standards Code, 2019 (CBSC 2019) for component supports and attachments are not permitted. For distribution system, interior partition wall, and suspended ceiling seismic bracings, test criteria other than those adopted in the CBSC 2019 may be used when approved by OSHPD prior to testing.							
X Analysis BY: Jeffrey Kikumoto							
Experience Data DATE: 05/06/2020							
Combination of Testing, Analysis, and/or Experience Data (Please Specify):							
COTÉ COTÉ DE COTE DE C							
OSHPD Approval BUILDING							
Date: 5/6/2020							
Name: Jeffrey Kikumoto Title: Senior Structural Engineer							
Condition of Approval (if applicable):							

"Access to Safe. Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"









5877 Pine Ave, Ste. 210 Chino Hills, CA. 91709 Phn: (909) 606-7622

Office of Statewide Health Planning and Development

PREAPPROVAL OF MANUFACTURER'S CERTIFICATION OPM-0208

THIS PREAPPROVAL CONFORMS TO THE 2019 CALIFORNIA BUILDING CODE

MANUFACTURER:

EQUIPMENT NAME:

BD

PYXIS DUOSTATION

Sheet: <u>1 of 12</u>


Date: 4/29/20

GENERAL NOTES

- 1. THIS OSHPD PREAPPROVAL OF MANUFACTURER'S CERTIFICATION (OPM) IS BASED ON THE 2019 CBC. THE DEMANDS (DESIGN FORCES) FOR USE WITH THIS OPM SHALL BE BASED ON THE 2019 CBC
- 2. THIS DOCUMENT MAY ONLY BE USED WITH THE EXPRESS WRITTEN CONSENT OF THE MANUFACTURER LISTED ABOVE FOR THE SPECIFIC PROJECT SITE AND INSTALLATION LOCATION. THIS DOCUMENT IS INVALID WITHOUT SUCH CONSENT.
- 3. THIS PREAPPROVAL CONFORMS TO THE 2019 CALIFORNIA BUILDING CODE WHERE SDS IS NOT GREATER THAN 1.00 & 1.85. SEE DETAIL FOR APPLICABILITY
- 4. FORCES PER ASCE 7-10 SECTION 13.3.1, EQUATIONS 13.3-1, 13.3-2 & 13.3-3,
 - WHERE SDS = 1.00, \mathbf{a}_{p} = 1.0, \mathbf{I}_{p} = 1.5, \mathbf{R}_{p} = 1.5, \mathbf{z}/h = 0 AT CONCRETE SLAB, $\mathbf{z}/h < 1$ AT CONCRETE SLAB ON METAL DECK SEE FOLLOWING SHEETS FOR Ω_{0}
 - WHERE SDS = 1.85, \mathbf{a}_p = 1.0, \mathbf{I}_p = 1.5, \mathbf{R}_p = 1.5, \mathbf{z}/h = 0.4T CONCRETE SLAB, $\mathbf{z}/h \le 1$ AT CONCRETE SLAB ON METAL DECK. SEE FOLLOWING SHEETS FOR Ω_0
- THIS PREAPPROVAL COVERS ONLY THE SUPPORTS AND ATTACHMENTS OF THE EQUIPMENT TO THE STRUCTURE.
- 6. ALL DESIGN FORCES SHOWN ON THE DRAWINGS ARE FACTORED LOADS THAT SHALL BE USED FOR STRENGTH DESIGN.
- 7. CONCRETE SLAB ON METAL DECK DETAIL VALID FOR DEMANDS SHOWN AT ANY ELEVATION IN THE BUILDING. (i.e. z/h < 1)
- 8. CONCRETE SLAB DETAIL VALID FOR DEMANDS SHOWN AT ANY ELEVATION AT OR BELOW GRADE. (i.e. z/h = 0)

9. RESPONSIBILITIES OF THE STRUCTURAL ENGINEER OF RECORD OF THE BUILDING

- A. PROVIDE SUPPORTING STRUCTURE TO SUPPORT WEIGHTS AND FORCES SHOWN IN ADDITION TO ALL OTHER LOADS.
- B. VERIFY THAT THE INSTALLATION IS IN CONFORMANCE WITH THE 2019 CBC AND WITH THE DETAILS, MATERIAL AND GAGE OF THE UNIT WHERE ATTACHMENTS ARE MADE AGREE WITH THE INFORMATION SHOWN ON THE PREAPPROVAL DOCUMENTS.
- C. VERIFY THAT PROJECT SPECIFIC VALUES OF SDS & z/h RESULT IN SEISMIC FORCES (Eh, Ev) THAT DO NOT EXCEED THE VALUES ON THE DETAILS.
- D. VERIFY THAT THE CONCRETE SLAB TO WHICH THE EQUIPMENT IS ANCHORED MEETS THI REQUIREMENTS OF THE APPLICABLE ICC ESR REPORT AND THIS OPM.
- E. VERIFY THAT THE ANCHORS ARE AN ADEQUATE DISTANCE FROM ANY SLAB EDGES OR OPENINGS (SEE TYPICAL DETAIL ON SHEET 2).
- F. VERIFY THAT ALL NEW OR EXISTING ANCHORS ARE AN ADEQUATE DISTANCE FROM THE UNIT ATTACHMENTS AND CHECK FOR INTERACTION WHERE OTHER ANCHORS ARE WITHIN 18" OR 6hef FROM THIS UNIT'S ANCHORS.



www.EquipmentAnchorage.com

DES. J. ROBERSON

DATE

јов no. 11**-2009**

4/29/20

. 12 «UEET»

SHEET

PYXIS DUOSTATION

BD

10. EXPANSION ANCHORS:


A. ATTACHMENT IS TO BE MADE WITH THE ANCHORS LISTED BELOW AND INSTALLED AS DESCRIBED IN THE CORRESPONDING ICC REPORT.

Anchor Diameter	Concrete Type	Min. f'c (psi)	Anchor Type	ICC Report No.	Min. Embed.	Min. Spacing	Min. Edge Dist	Min. Conc. Thickness	Torque Test	Direct Tension
1/2"	Sand Light Weight	3000	Hilti Kwik Bolt TZ	ESR-1917	3.25"	6.75"	24"	See Detail "A"	40 FT-LB	1585 lb
1/2"	Normal Weight	3000	Hilti Kwik Bolt TZ	ESR-1917	3-1/4"	8"	24"	6"	40 FT-LB	3026 lb


- B. THIS PREAPPROVAL ALLOWS FOR UP TO A MAXIMUM OF 2 ADJACENT CONCRETE SLAB EDGES, 24" AWAY MINIMUM (i.e. CORNER). SEE ADJACENT DETAIL FOR ADDITIONAL MINIMUM ALLOWABLE CONCRETE EDGE DISTANCES.
- C. TESTING OF EXPANSION ANCHORS PER 2019 CBC, 1913A.7:
 TESTING SHALL BE DONE IN THE PRESENCE OF THE SPECIAL
 INSPECTOR AND A REPORT OF THE TEST RESULTS SHALL BE
 SUBMITTED TO OSHPD
 - (i) AFTER AT LEAST 24 HOURS HAVE ELAPSED SINCE INSTALLATION,
 DIRECT PULL TENSION TEST OR TORQUE TEST AT LEAST 50% OF THE ANCHORS.
 - (ii) ACCEPTANCE CRITERIA:

DATE: 05/06/2020

- DIRECT TENSION TEST: THE ANCHOR SHOULD HAVE NO OBSERVABLE MOVEMENT AT THE TEST LOAD. A PRACTICAL WAY TO DETERMINE OBSERVABLE MOVEMENT IS THAT THE WASHER BECOMES LOOSE.
- TORQUE TEST: THE APPLICABLE TORQUE MUST BE ACHIEVED
 WITHIN THE FOLLOWING LIMITS: WEDGE TYPE: 1/2 TURN OF THE
 NUT
- (iii) IF ANY ANCHOR FAILS. TEST ALL ANCHORS.
- 11. BOLTS THROUGH CONCRETE ON METAL DECK
 - A. BOLTS SHALL BE TORQUED BY 3/4 TURN OF THE NUTS AFTER THE SNUG TIGHT (THE SNUG-TIGHT CONDITION IS DEFINED AS THE TIGHTNESS REQUIRED TO BRING THE CONNECTED PLIES INTO FIRM CONTACT) CONDITION IS ACHIEVED, UNLESS OTHERWISE NOTED,
 - B. THROUGH BOLT HOLES SHALL BE 1/16" LARGER THAN BOLT SIZE (HOLE SIZE = BOLT SIZE + 1/16) FOR CONCRETE.
 - C. THROUGH-BOLTS IN CONCRETE SHALL RECEIVE SPECIAL INSPECTION AND TESTING (THROUGH BOLTS WITH STEEL TO STEEL CONNECTION IN TENSION DO NOT REQUIRE TENSION TESTING) IN ACCORDANCE WITH REQUIREMENTS FOR POST-INSTALLED ANCHORS.



TYPICAL CONCRETE EDGE DETAIL



www.EquipmentAnchorage.com

BD

PYXIS DUOSTATION

DES. J. ROBERSON


JOB NO. 11-2009

3

DATE 4/29/20

of 12 SHEETS

SHEET



NOTES:

1. FORCES ARE DETERMINED PER 2019 CALIFORNIA BUILDING CODE AND ASCE 7-16.

STRENGTH DESIGN IS USED. (ap = 1.0, lp = 1.5, Rp = 1.5, Ω_0 = 1.5, z/h = 0)

2. CENTER OF GRAVITY (C.G.) AND WEIGHT ARE THE GOVERNING PARAMETERS FOR DESIGN. THIS PREAPPROVAL ENCOMPASSES ALL WEIGHTS UP TO THE MAXIMUM WEIGHT SHOWN.

3. STRUCTURAL ENGINEER OF RECORD FOR THE BUILDING SHALL PROVIDE SUPPORT STRUCTURE DESIGNED TO SUPPORT WEIGHTS AND FORCES SHOWN IN COMBINATION WITH ALL OTHER LOADS THAT MAY BE PRESENT.

4. CALCULATION COVERS PYXIX STATION DOUBLE UNIT, SYSTEM 30 DOUBLE UNIT. & SYSTEMS EC DOUBLE UNIT. ALL UNITS HAVE SAME WEIGHT AND DIMENSIONS.



EASE

EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

www.EquipmentAnchorage.com

BD

PYXIS DUOSTATION

DES. J. ROBERSON


JOB NO. 11-2009

4

SHEET

DATE 4/29/20

OF 12 SHEETS





EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING www.EquipmentAnchorage.com DES. J. ROBERSON BD

PYXIS DUOSTATION


11-2009 JOB NO.

DATE

4/29/20

SHEET

SHEETS





www.EquipmentAnchorage.com

BD

DES. J. ROBERSON

DATE

11-2009 JOB NO.

4/29/20

SHEET

SHEETS

PYXIS DUOSTATION

28.3" FOR SINGLE SEISMIC SUPPORTS & ATTACHMENTS CONCRETE SLAB ON METAL DECK 51.75" FOR DOUBLE C.G.WT. = (SEE SCHEDULE) (SEE SHEET 9 OF 12) SEISMIC ANCHORING KIT (PROVIDED BY CAREFUSION) SEE DETAILS "B" AND "C" USE 4- 1/2" (A449) THREADED RODS THRU FLOOR N.W. OR SAND L.W. CONC. (3000 PSI MIN.) REFER TO MIN STEEL DECK REQUIREMENTS AND STRUT DETAIL "A" FRONT ELEVATION FOR ADDITIONAL DETAILS BUILDING

NOTE: SEE SHEET 7-8 OF 12 FOR WEIGHTS, DIMENSION & BOLT FORCES

NOTES:

FORCES ARE DETERMINED PER 2019 CALIFORNIA BUILDING CODE AND ASCE 7-16

STRENGTH DESIGN IS USED. (ap = 1.0, lp = 1.5, Rp = 1.5, Ω_0 = 1.5 z/h < 1)

Sds	1.00	1.85
HORIZONTAL FORCE (Eh)	1.20 Wp	2.22 Wp
HORIZONTAL FORCE (Emh)	1.80 Wp	3.33 Wp
VERTICAL FORCE (Ev)	0.20 Wp	0.37 Wp

 $(Emh = Eh \times \Omega_0; FOR CONCRETE ANCHORAGE)$

2. CENTER OF GRAVITY (C.G.) AND WEIGHT ARE THE GOVERNING PARAMETERS FOR DESIGN, THIS PREAPPROVAL ENCOMPASSES ALL WEIGHTS UP TO THE MAXIMUM WEIGHT SHOWN.

3. STRUCTURAL ENGINEER OF RECORD FOR THE BUILDING SHALL PROVIDE SUPPORT STRUCTURE DESIGNED TO SUPPORT WEIGHTS AND FORCES SHOWN IN COMBINATION WITH ALL OTHER LOADS THAT MAY BE PRESENT.



EASE

EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

www.EquipmentAnchorage.com

BD

DES. J. ROBERSON

JOB NO. 11-2009


7

SHEET

PYXIS DUOSTATION

DATE 4/29/20

_F 12 _{энеетэ}





BD

EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

www.EquipmentAnchorage.com

DES. J. ROBERSON

11-2009 JOB NO.


SHEETS

SHEET

PYXIS DUOSTATION

4/29/20 DATE

SEISMIC SUPPORTS & ATTACHMENTS MAX Sps ≤ 1.85 CONCRETE DETAIL USE 4- 1/2" (A449)



PLAN AT BASE (SINGLE 3-DRAWER) BUILDING

MODEL Eı (in) E2 (in) SINGLE 3 DWR 12.83 13.67 SINGLE 6 DWR 12.40 14.10



(EA SIDE)

www.EquipmentAnchorage.com

BD

DES. J. ROBERSON

11-2009

SHEET

PYXIS DUOSTATION

4/29/20 DATE

JOB NO.

OF SHEETS

SEISMIC SUPPORTS & ATTACHMENTS

BASE DETAIL

	SDS	MODEL	WEIGHT (LB)	₹ (in)	*** T U (lb)	*** (lb)
**	1.85	SINGLE 3-DWR	1425	42.94	3857	953
	1.25	SINGLE 6-DWR	1338	45.1	3814	895
*	1.00	DOUBLE 3-DWR	R28730	45.07	4100	III8
	1.00	DOUBLE 6-DWR	2810	45.59	4063	1093

* THIS UNIT IS USED IN CALCULATION ON SHEET 5 OF 10 ** THIS UNIT IS USED IN CALCULATION ON SHEET 6 OF 10

*** (VALUES DO NOT INCLUDE 10.) 208

No. 4197

EASE

EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

www.EquipmentAnchorage.com

BD

DE8. J. ROBERSON

JOB NO. 11-2009


10


SHEET

PYXIS DUOSTATION

DATE 4/29/20

of 12 SHEETS





BD PYXIS DUOSTATION

EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

WWW. EquipmentAnchorage.com

SHEET

DES. J. ROBERSON

JOB NO. 11-2009

4/29/20

DATE

OF 12 SHEETS


SEISMIC SUPPORTS & ATTACHMENTS

RIGHT & LEFT
FLOOR BRACKETS
(I/4" THK, ASTM A512,
GRADE 50, Fy=50 KSI (MIN))

0.75" 8.25" | " 0.50" 26.50"

<u>SECTION A-A</u>

SECTION B-B



BRACKET DETAIL

www.EquipmentAnchorage.com

BD

DES. J. ROBERSON

4/29/20

11-2009 JOB NO.

DATE

SHEETS

SHEET

PYXIS DUOSTATION

SEISMIC SUPPORTS & ATTACHMENTS BASE DETAIL STEEL DOWEL (1/2" DIA, ASTM AIOI8 HSLAS, CLASS I, GRADE 50, Fy=50 KSI (MIN)) (EA SIDE) Ζ "a" "b" STEEL DOWEL 0 C.G. Χ OPM-0208 1/2" STEEL PLATE (ASTM AIOI8 HSLAS, CLASS I <u>_</u>o GRADE 50, Fy=50 KSI (MIN)) Jeffrey Kikumoto FRONT FACE OF BRACKET

(4)- 1/2"Φ CAP SCREWS (GRADE 5) TO TAPPED HOLES AT UNIT (SEE SCHED)

NOTE: THE SEISMIC ANCHORING KIT CONSISTS OF TWO FLOOR BRACKETS, A RIGHT FLOOR BRACKET AND A LEFT FLOOR BRACKET (SEE DETAIL "B") AND FRAME ASSEMBLY (SHOWN ABOVE) (TYP)

REFLECTED PLAN AT FRAME ASSEMBLY DETAIL

(FRAME ASSEMBLY SHOP ATTACHED TO UNDERSIDE OF UNIT BY MFG)

MODEL	"a" (in)	"b" (in)	"c" (in)	"d" (in)	"e" (in)	"f" (in)
DOUBLE 6-DWR	22.90	25.82	22.46	25.38	10.95	12.03
DOUBLE 3-DWR	22.90	25.82	22.46	25.38	10.95	12.03
SINGLE 6-DWR	12.35	14.32	11.91	13.88	12.36	10.62
SINGLE 3-DWR	12.35	14.32	11.91	13.88	12.36	10.62

